

LipiORDER™ <Membrane Lipid Order Imaging Dye>

Catalog NO. FDV-0041

Research use only, not for human or animal therapeutic or diagnostic use.

日本語版はこちらから
ダウンロードできます。
①弊社ウェブサイトより
Webページ番号検索にて
【70873】で検索



Product Background

Membrane lipid order is a biophysical parameter that defines a membrane organization and is often described by the degree of lipid packing. For example, phospholipids only containing saturated lipids create high packing and thick lipid bilayer, called liquid-order (Lo) phase. On the other hand, phospholipids containing unsaturated lipids, which have bent structure, form low packing and thin membrane structure, called liquid-disorder (Ld) phase. In the model membrane mixing saturated lipids and unsaturated lipids, Lo and Ld are clearly separated and create individual domains. While the model membrane composition can be discussed membrane lipid order (Lo/Ld) easily, actual cells have numerous types of lipids and form very complicated membrane lipid orders. Furthermore, the lipid order is also influenced by various factors in cellular membranes, including sterol lipids such as cholesterol and membrane proteins, etc. Lipid raft, a continuous interest topic in biology, which serves as functional microdomains on cellular membranes, is one of the specialized Lo domains, with highly accumulated saturated lipids such as sphingomyelin, cholesterol, functional membrane proteins and lipidated proteins. Membrane lipid order has been considered as a fundamental factor in providing physical properties of cellular membranes, such as membrane fluidity, membrane tension and membrane curvature. Observation of cellular lipid order may lead to an understanding of the various function of cellular membranes.

To measure membrane lipid order, some solvatochromic dyes which change fluorescence intensity and color in response to their solvent polarity are applied. These solvatochromic dye fluorescent properties change depending on membrane lipid order. Among them, Laurdan is the most well-known dye for membrane lipid order imaging. However, conventional dyes have some limitations. For example, Laurdan requires UV light excitation and exhibits low photostability. So Laurdan is not suitable for live-cell imaging. Dyes which can be excited by longer wavelength with more photostability and chemically stable in cells are desirable traits for cellular imaging of membrane lipid order. LipiORDER™ is a novel solvatochromic dye for membrane lipid order imaging (original compound name PK in Ref.1). LipiORDER™ is excited at around 400 nm wavelength, which is compatible with live-cell imaging and changes its emission fluorescent color from green to red depending on membrane lipid order. LipiORDER™ also has high photostability and chemical stability on the cell membranes. LipiORDER™ is a convincing tool to monitor cellular membrane lipid order imaging on live-cells.




Figure 1. Overview of membrane lipid order

Principle and Reference data

Sensing of lipid order by using LipiORDER™ is based on the following two unique properties. 1) LipiORDER™ is a pyren-based solvatochromic fluorescent dye which changes fluorescent property in response to their solvent environment (Figure P1). In low polaric solvents such as toluene, LipiORDER™ shows green fluorescence. On the other hand, in highly polaric solvents such as DMSO and methanol, this dye changes color to orange or red. 2) LipiORDER™ is a highly hydrophobic compound and quickly accumulates in the various biological membranes. Combining the two features above, LipiORDER™ can sense the local environment in a lipid bilayer. Generally, Lo is a high packing lipid bilayer and shows lower polarity, whereas Ld is a sparse packing lipid bilayer and shows high polarity. Based on polarity of lipid bilayer derived from lipid order, LipiORDER™ will change fluorescent color, from green on Lo membrane to red on Ld membrane (Figure P2). Ratiometric fluorescent value (F_{Red}/F_{Green}) is correlated to lipid order (Lo and Ld).

Actually, in sphingomyeline/cholesterol (SM/Chol) liposome, one of the model Lo, LipiORDER™ emits green fluorescence and in 1,2-dioleoyl-sn-glycerol-3-phosphocholine (DOPC) liposome, a model Ld, shows red fluorescence. In DOPC/Chol, an intermediate model, the reagent show yellow to orange. The ratiometric values (F_{575}/F_{510}) clearly depend on lipid order, SM/Chol (Lo) is low and DOPC (Ld) is high (Figure P3).




Figure P1 Absorption and fluorescent spectrum of LipiORDER™ in various solvent




Figure P2 Graphical overview of lipid order-dependent fluorescent change of LipiORDER™




Figure P3 Fluorescent spectrum of LipiORDER™ in model liposomes

Description

Catalog Number: FDV-0041

Size: 0.1 mg

Formulation: C₂₃H₂₁NO

Molecular weight: 327.4 g/mol

Solubility: Soluble in DMSO

Fluorescent characteristics: Ex. 405 nm/Em. 450-650 nm (dependent on solvents)

Reconstitution and Storage

Reconstitution: Stock solution recommended concentration 1 mM in 100% DMSO.

Storage (powder): Store powder at -20°C.

Storage (solution): After reconstitution in DMSO, aliquot and store at -20 °C.

Avoid repeated freeze-thaw cycles.


How to use and experimental setting

General procedure of live cell imaging

*This procedure is an example of cultured cell staining. For zebrafish staining, please find Ref.1.

1. Prepare 0.1-1 μM LipiORDER™ in serum-free and phenol red-free medium such as HBSS
2. Remove culture medium and wash cells PBS several times
3. Add LipiORDER-containing medium to cells
4. Incubate cells at 37 °C for over 10 min
5. Wash cells with PBS or medium (Optional)
6. Observe cells under live condition with confocal laser microscopy and obtain green and red fluorescent images
7. Perform ratiometric analysis with image software using green and red fluorescent images

NOTE: The staining concentration of LipiORDER™ is dependent on cell type and experiments. Please empirically optimize to determine the suitable concentration for each experiment.



Fluorescent microscopy and analysis

For LipiORDER™ ratiometric imaging, LipiORDER™ is excited at 405 nm and its fluorescence is detected with two ranges, green channel and red channel. The recommended wavelength range of green channel and red channel is 500-550 nm and 550-650 nm, respectively. Ratiometric image analysis (F_{Red}/F_{Green}) is calculated by any image processing software such as ImageJ.

Option: For calibration control of each model lipid order, we recommend obtaining the following images.

- Labrafac oil for lipid droplet,
- SM/Chol liposome for Lo model
- DOPC liposome for Ld model


NOTE: As you can see in Figure P1, the absorbance of LipiORDER™ at around 480 nm is negligible.

LipiORDER™ is compatible with common green dyes (excited by ~480 nm laser) and red dyes (excited by ~560 nm) for multicolor staining.

Application data


Photostability of LipiORDER™

LipiORDER™ and Laurdan, a conventional membrane lipid order imaging dye in lipid vesicles composed of 0.2 mM DOPC in 20 mM HEPES (pH 7.4) were irradiated with Xe lamp. LipiORDER™ and Laurdan were excited at 405 nm and 360 nm, respectively and fluorescent intensity was measured. Laurdan was quickly photodegraded, whereas LipiORDER™ maintains fluorescent intensity for at least 1 hour. LipiORDER™ is highly stable compared to Laurdan.




Ratiometric imaging of COS7 cells

COS7 cells were treated with 300 nM LipiORDER™ in HBSS for 10 min and observed by confocal laser microscopy (Ex. 405 nm, Em 470-550 nm for Green channel and >550 nm for Red channel). Ratiometric analysis was performed with ImageJ using green and red channel data and lipid order was shown by green-to-red pseudocolor (Lo Ld). Plasma membrane and intramembranes are shown Lo and Ld, respectively.




*An example of membrane lipid order analysis

Each layer of ratiometric pseudocolor was extracted. Low ratio value shows plasma membrane structure mainly and high ratio value and mainly shows how ER-like and mitochondria-like structure, respectively.




Ratiometric imaging of neuronal cells


Primary cultured hippocampal neurons (DIV 3 or DIV 12) from E17.5 mice were stained with 300 nM LipiORDER™ in HBSS for 10 min and observed by confocal laser microscopy (Ex. 405 nm, Em. 470-550 nm for Green channel and >550 nm for Red channel). Ratiometric analysis was performed with ImageJ using green and red channel data and lipid order was shown by green-to-red pseudocolor (Lo Ld).





* Each fluorescent imaging data is shown below. The ratiometric data was calculated using the following pictures.




*An example of membrane lipid order analysis (Dendrites in DIV12 neurons). Near Lo phases () clearly shows plasma membrane structures. On the other hand, intermediate () and Ld phases () were observed from intracellular compartments.



Drug-induced cellular lipid order changes

COS7 cells were treated with 15 mM β -cyclodextrin (β -CD), a membrane-disrupting chemical via removing endogenous cholesterol, for 4 hours. After β -CD treatment, cells were washed and stained with 300 nM LipiORDER™ in HBSS for 10 min. The cells were observed by confocal laser microscopy (Ex. 405 nm, Em. 470-550 nm for Green channel and >550 nm for Red channel). Ratiometric analysis was performed with ImageJ using green and red channel data and lipid order is shown by green-to-red pseudocolor (Lo  Ld). The cell structure was dramatically changed by β -CD and at the same time, the distribution of Lo phase () clearly changed.



Notes

All spectrum data and a photostability data were obtained by Dr. Yosuke Niko, Kochi University.

All cell imaging data were obtained by Dr. Mitsuharu Hattori, Nagoya City University.

Reference

1. Valanciunaite *et al.*, *Anal. Chem.*, **92**, 6512-6520 (2020) Polarity Mapping of Cells and Embryos by Improved Fluorescent Solvatochromic Pyrene Probe.

Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter “Product information”) have been prepared based on published research reports on November, 2020. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本商品は学術研究成果を基にフナコシ株式会社が商品化したもので、2020年11月時点における公開研究報告を基に広告文章およびマニュアル(以下、商品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、商品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本商品の仕様および商品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただけますようお願い申し上げます。

 E-mail Newsletter Sign Up

Japanese 

English 

 **funakoshi**
FRONTIERS IN LIFE SCIENCE

URL: <http://funakoshi.co.jp>
9-7 Hongo 2-Chome, Bunkyo-ku, Tokyo 113-0033

Related products

Lipidye™ II <Lipid Droplet Live Imaging>


Lipidye™ II is a highly sensitive lipid droplet staining dye with extremely photostable property. This dye is the second generation of our previous reagent, Lipidye™. This dye allows us to detect small lipid droplets (<1 μm) in non-adipocytes and to apply into long-term live cell imaging for dynamic lipid droplet movements.

Catalog No. FDV-0027

Size 0.1 mg

Features

- Recommended Ex/Em: 400-500 nm / 490-550 nm
- Enable to detect <1 μm lipid droplets
- Suitable for long-term live cell imaging
- Extremely photostable compared with conventional dyes
- Compatible with both live and fixed cells



FAOBlue™ <Fatty Acid Oxidation Detection Reagent>


FAOBlue™ is a cell-based fatty acid beta-oxidation (FAO) detection dye which emits blue fluorescence upon FAO activity. FAOBlue™ enables to quantitatively monitor cellular FAO activities under various conditions.

Catalog No. FDV-0033

Size 0.2 mg

Features

- Recommended Ex/Em: ~405 nm / 460 nm
- Enable to detect cellular FAO activity directly without any specific equipment, only need microscopy.
- Monitor drug-induced change of FAO activity quantitatively.



LipirADICAL™ Green <Lipid Radical Detection Reagent>

LipirADICAL™ Green is a specific fluorescent dye for lipid-derived radicals which are the most upstream factor of lipid peroxidation (LPO). LipirADICAL™ Green can be applied into both *in vitro* assay and cell-based assay to monitor lipid radical productions.

Catalog No. FDV-0042

Size 0.1 mg

Features

- Recommended Ex/Em: ~480 nm / 520 nm
- Enable to detect very unstable lipid-derived radicals
- Compatible with *in vitro* assay and in cell-based assay
- An innovative reagent for comprehensive identification of lipid-derived radicals by lipidomics

