


ZF Plate Instruction Manual

The ZF plate is a 96-well microplate with a unique well-slit design that allows efficient alignment of zebrafish. This microplate with optical clear bottom enables high-throughput imaging for zebrafish.

1. Package contents

- ZF Plate (with transparent cover): 5 plates.
- Instruction Manual: One

2. Component descriptions



3. Procedure

1. Required materials and equipment:
 - Anesthetized zebrafish (approximately 3dpf to 6dpf)
 - ZF plate
 - Pipets (100 ul to 200 ul and 1 ml)
 - Pipet tips (wide bore tip recommended)

- 1) Anaesthetized zebrafish are placed into the wells
 - Using a wide-bore tip, gently pipet a single anesthetized zebrafish and place into an available well, repeat as necessary.
- 2) Quick centrifuge
 - Using the flash function, quickly centrifuge the microplate at 200 to 250Gs for two seconds to place the zebrafish within the slit.
- 3) Image zebrafish
 - After imaging zebrafish, add additional anesthetic water into the wells and using a pipet with a wide-bore tip gently suction out the zebrafish from each well.

See the URL below for the more details regarding the procedure.

<https://www.hashimoto-inc.co.jp/products/zfプレート/>



4. Specification

Model	HDK-ZFA101-02a	Plate material/color	polystyrene / black
Plate size(mm)	85.48(W)×127.76(L)×14.35(H)	Bottom material/thickness(mm)	Silica glass / 0.17
Mass(g)	without cover : 72 / with cover : 95	Temperature for use	0 ~ 40°C
Applicable standards	ANSI/SLAS-1 ~ 4 (96well)	Temperature for storage	4 ~ 25°C
Per well maximum capacity(ul)	250	Sterilized/non-sterilized	non-sterilized

5. Precautions for use

- ⚠ Do not use for applications other than imaging of small fish.
- ⚠ This component is disposable. Do not reuse.
- ⚠ Cultivation of small fish in the well is not possible.
- ⚠ This product cannot be treated with an autoclave.



HASHIMOTO ELECTRONIC INDUSTRY CO.,LTD
3866-12, Takasucho, Matsusaka-shi, Mie, 515-0104, Japan
TEL:+81(0)598-51-3111 FAX : +81(0)598-52-1417
<https://www.hashimoto-inc.co.jp/>