

# iPSelector < Anti-LNFP I, Human, Mouse-Mono (R-17F)>

Catalog NO. FDV-0014B

Research use only, not for human or animal therapeutic or diagnostic use. This product has been commercialized under the license from Ritsumeikan University.

## **Product Background**

Clone SSEA-3, SSEA-4, TRA-1-60, and TRA-1-81 antibodies are well-known as human iPS/ES cell-marker antibodies. Since SSEA-3 antibody was originally generated against mouse embryo and the other antibodies were against human EC cells, these antibodies recognize not only human iPS/ES cells but also human EC cells.

Our **iPSelector** (**clone R-17F**) is a novel mouse monoclonal antibody generated by using a human iPS cell line as an immunogen. It is specific to human iPS/ES cells and does not essentially cross-react against human EC cells (Table 1, ref. 1). This iPSelector (clone R-17F) antibody also stains entire surface of human iPS/ES cell membranes evenly, while the staining by SSEA-3 and SSEA-4 antibodies are not uniformly (ref. 2). In addition, iPSelector (clone R-17F) is reported to exhibit potent dose-dependent cytotoxicity when added to living human iPS/ES cells (ref. 2 & 3). iPSelector (clone R-17F) is a beneficial tool for the selective detection, staining and removal of undifferentiated of human iPS/ES cells in regenerative medicine.

Table 1. Binding Activity of Antibodies to Cells

	iPSlector R-17F	TRA-1-60	TRA-1-81	SSEA-3	SSEA-4
Tic (iPS)	++++	++++	++++	++++	++++
KhES-3 (ES)	+++	++++	++++	+++	++++
H9 (ES)	++++	++++	++++	+++	++++
2102Ep (EC)	+/-	++++	++++	+++	+++

# **Description**

Catalog Number: FDV-0014B

Size: 100 µL

Lot No.: see vial label

Host Species and Clonality: Mouse Monoclonal

Clone name: R-17F

Specificity: This antibody recognizes lacto-N-fucopentaose I (LNFP I: Fucα1–2Galβ1–3GlcNAcβ1–3Galβ1–

4Glc) on a glycolipid / glycoprotein. R-17F epitopes are expressed on undifferentiated human induced pluripotent stem (iPS) / embryonic stem (ES) cells but not on human embryonal carcinoma

(EC) cells nor on differentiated human iPS/ES cells.

Isotype and Subclass: IgG1

Formulation: Phosphate Buffered Saline (PBS) containing 50% Glycerol, contains no preservative.

Purification: Protein A Purified Concentration: 1.0 mg/ml

Verified Species Reactivity: Human \* Note: Other species not tested.

Immunogen: Human iPS cell line, Tic, derived from human fetus lung cells (MRC-5).

Storage: -20°C (Avoid repeated freeze-thaw cycles.)


### **Application and Recommended usage**

- Western blotting 1/2,000

Immunohistochemistry
Flow cytometry
Functional applications
Optimal dilutions should be empirically determined for each experiments
Optimal conditions should be empirically determined for each experiments

# **Reference and Application Data**

#### Western Blotting



Sample: 5 µg cell lysate in each lane. Left: human iPS cells (LNFP I positive) Right: HEK293 (Negative Control)

**Dilution:** 1:2,000

Secondary antibody: Anti-Mouse IgG, Goat-Poly, HRP

(Kirkegaard & Perry Laboratories, #5220-0337)

Chemiluminescence Substrate: Trident plus Western HRP Substrate


(GeneTex, #GTX400006)

Detection: LuminoGraph I (ATTO) with 1min exposure

One major positive band and several minor bands were specific to human iPS cells, and any positive band was not obtained with HEK293 cells.

#### Immunocytochemical Staining

Cultured human iPS cells were stained with R-17F, SSEA-3, and SSEA-4 antibodies. [bars:  $100 \mu m$ ] R-17F stained the entire surface of the cell membranes equally, while the staining by SSEA-3 and SSEA-4 antibodies are not evenly. This suggests that R-17F epitope is expressed ubiquitously all over the human iPS cells




#### Functional Application (Cytotoxic effects on undifferentiated iPS/ES cells)

R-17F is reported to exhibit potent dose-dependent cytotoxicity when it is added to living human iPS/ES cells.

[Left] After the incubation of iPS cell suspension with R-17F at 4°C for only 45 minutes, the percentage of viable cells decreased concentration-dependently (red triangles).

Blue squares: effects of the isotype (IgG1)-matching control antibody (anti- $\alpha$ -MBP) as Negative Control [Right] When R-17F-treated iPS cells were incubated with a small amount (0.025-0.1  $\mu$ g) of the secondary antibody (goat anti-mouse IgG1 antibody), the cytotoxic effect of R-17F was enhanced significantly in a dose-dependent manner (red triangles).



#### Reference

- 1. Kawabe, *et al.*, *Glycobiology*, **23**, 322-336 (2013) A novel antibody for human induced pluripotent stem cells and embryonic stem cells recognizes a type of keratan sulfate lacking oversulfated structures
- 2. Matsumoto, *et al.*, *J. Biol. Chem.*, **290**, 2007-200851 (2015) A Cytotoxic Antibody Recognizing Lacto-N-fucopentaose I (LNFP I) on Human Induced Pluripotent Stem (hiPS) Cells
- 3. Nakao, *et al.*, *Glycoconj. J.*, **34**, 779-787 (2017) Characterization of glycoproteins expressing the blood group H type 1 epitope on human induced pluripotent stem (hiPS) cells



