

Introduction

Trajan emitters are optimized for high resolution with extended life in a wide variety of mass spectrometry instrumentation, applications, and workflows.

Intended use

Trajan's nanospray emitters (ID: 10 - 20 μm) are designed specifically for low-flow applications and offer premium performance and robustness/longevity.


Important

⚠ Precautionary statements:

- Caution should be taken to avoid injury (especially to the eyes) when working with fused silica (quartz). It is recommended to use safety glasses and avoid unnecessary bending or breaking of the fused silica.
- In the event of an emitter breakage or emitter replacement, observe local requirements for disposal of 'broken glass'.
- Narrow-bore emitters can lead to higher back-pressure in the system if partial blockages occur; be careful not to over-pressurize the system as the emitter may dislodge from the fitting/connection.
- Overtightening of the fused silica emitter can cause fractures, leading to higher risk of dead-volume and instrument contamination.

Instructions

1. It is recommended to avoid using metallic tweezers when handling emitters; this can cause irreversible damage to the fused-silica and/or metallic coatings (where applicable). expere voloritatent offic tem et qui cum dessimos soluptur, iumque. senimusa volut architis pa corum expedias nonsequid magniatis quo etus.



2. Open the lid of the emitter packaging, tearing the perforations in the label.



3. Using the appropriate PPE, carefully remove the desired emitter from the packaging, tip end first, using caution to avoid damaging the metal-coated and tip area of the emitter.



4. Using caution, install the emitter using the appropriate fittings, ferrules, and/or connections that correspond to your instrument configuration being careful not to overtighten and damage the fused-silica.

Distributed by:

CliniSciences Group Page 1 of 02

Storage and stability

Packaging should be stored at room temperature and ideally in a clean environment.

It is recommended to close the packaging lid when not actively retrieving an emitter to avoid contamination of the emitters with airborne dust and debris.

Troubleshooting

1. For the majority of issues, please first refer to the instrument troubleshooting guide.
2. For drops in performance that may be related to the emitter, try the following:
 - a) Turn off the voltage, check for a droplet accumulating at the tip. If yes, the flow is not blocked. Discard if blocked.
 - b) If there is no droplet at the tip with the voltage
3. If emitter issues are unable to be resolved, please contact tech support on techsupport@trajanscimed.com.

Disposal

Refer to local regulations for disposal in an approved waste disposal facility or service. Please recycle the packaging where possible.

Information and support

Refer to local regulations for disposal in an approved waste disposal facility or service. Please recycle the packaging where possible.

Distributed by:

CliniSciences Group

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: danmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre - France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

