


Stock Name : Novoprotein

Stock Code: 688137

**novoprotein**  
近岸蛋白

# GMP Grade Products



Novoprotein Scientific Inc.

**Distributed by:**

**CliniSciences Group**

# INDEX

<b>1</b>	Novoprotein's GMP Quality Management System	2
<b>2</b>	Regulatory Basis	2
<b>3</b>	Quality Control	2
<b>4</b>	FDA DMF	3
<b>5</b>	Files Support for IND	3
<b>6</b>	GMP Grade Products	4
	(1) Recombinant anti-Human CD3 mAb/CD28 mAb	4
	(2) NovoNectin®	5
	(3) Vitronectin	6
	(4) Cell/Organoid Culture Cytokines and Antibodies	7
	(5) Recombinant Cas9 Nuclease	9
	(6) (High Salt Active) BenzoNuclease®	10
	(7) mRNA Raw Materials and Reagents	14
<b>7</b>	GMP Grade Products Summary	17

Distributed by:

**CliniSciences Group**



## Dedicated & Professional


Novoprotein Scientific Inc. (Novoprotein) is a high-tech enterprise with more than 10 years of extensive experience in the recombinant protein industry, focusing on protein technology, and advanced in R&D, production, sales, and application solutions to raw materials and techniques for biopharmaceuticals, in vitro diagnosis, mRNA vaccines, and basic life science research. Our principal products include target proteins and cytokines, recombinant antibodies, molecular enzymes and reagents, as well as providing related technical services. Novoprotein possesses R&D and manufacturing bases in Shanghai, Suzhou, and Heze.

**novoprotein**  
近岸蛋白

**Distributed by:**

**CliniSciences Group**

# GMP Grade Key Raw Materials



**GMP Grade Factory** 8000m<sup>2</sup>

**FDA DMF** 15

Raw Materials for mRNA Drugs: 10

Raw Materials for Biologic Drugs: 5

**Gene and Cell Therapy Customers** 200+

**Help Customers Succeed in CDE/FDA IND or EUA** 16

Gene and Cell Therapy Customers: 8

mRNA Vaccine Customers: 8

Recombinant Protein Vaccine Customer: 1

Distributed by:

**CliniSciences Group**

# 1 Novoprotein's GMP Quality Management System

The GMP products provided by Novoprotein have been certified by ISO 9001:2015 and are produced in accordance with GMP quality management system. The GMP quality management system is implemented to ensure the controllability and traceability of materials and processes.

- Adequate number of personnel and a well-established personnel training process
- The factory facilities and equipment comply with GMP and are regularly maintained and verified
- A comprehensive document management system, including process procedures, quality standards and inspection methods, complete batch production and inspection records, etc
- A comprehensive supplier management process, raw material and product management process
- Strict production process management, effective pollution prevention and cross contamination measures
- A comprehensive quality assurance and control system, with each product audited by QC and released by QA
- Fully validated equipment, processes, and inspection methods



## 2 Regulatory Basis

- (1) 《Good Manufacturing Practice of Medical Products (2010 version, NMPA)》
- (2) 《GMP Appendix – Cellular therapeutic product》 National Medical Products Administration
- (3) 《Pharmacopoeia of the People’s Republic of China》 2020 version, National Pharmacopoeia Commission
- (4) USP Chapter <1043>, Ancillary Materials for Cell, Gene, and Tissue-Engineered Products
- (5) USP Chapter <92>, Growth Factors and Cytokines Used in Cell Therapy Manufacturing
- (6) Ph. Eur. General Chapter 5.2.12, Raw Materials of Biological Origin for the Production of Cell-based and Gene Therapy Medicinal Products

## 3 Quality Control

- Animal-free materials and pharmaceutical-grade excipients are used for the production
- The production process has virus inactivation and removal steps, and some products can provide virus removal process validation reports
- Conduct comprehensive testing of production cell lines in accordance with the requirements of the pharmacopoeia, ensuring no virus or external factor contamination, and provide inspection reports
- Sterile, ampicillin-free, mycoplasma-free
- Medicinal penicillin bottle packaging

Distributed by:

**CliniSciences Group**

## 4 FDA DMF

Drug Master Files (DMFs) are archive files submitted by the holder to the FDA, which contain detailed information on the production facilities, process flow, quality control, raw materials, packaging materials, and other processes used in the production, processing, packaging, and storage of drug products for human use. DMF can be used for one or more clinical applications (IND), innovative drug applications (NAD), simplified new drug applications (ANDA), bioproduct license applications (BLA), as well as amendments and supplements to various applications mentioned above. The raw material supplier submits the required technical content directly to the FDA for filing and obtains a filing number in the form of a DMF document. Drug applicants can use the DMF filing number directly instead of providing specific information about raw materials and excipients during the filing process.

**Novoprotein has obtained FDA DMF registration for 15 GMP Grade products, including 10 mRNA vaccine drug raw materials and 5 biopharmaceutical raw materials.**

Cat. No.	Product Name	DMF No.
GMP-C013	Recombinant Human IL-2	037809
GMP-CD47	Recombinant Human IL-7	037735
GMP-C016	Recombinant Human IL-15	037736
GMP-1701	Recombinant Cas9 Nuclease, GMP Grade	039085
GMP-1707	BenzoNuclease®, GMP Grade	035864
GMP-RE026	BsaI, GMP Grade	037810
GMP-RE057	BspQI, GMP Grade	038671
GMP-RE015	XbaI, GMP Grade	038881
GMP-E121	T7 RNA Polymerase, GMP Grade	038343
GMP-E125	RNase Inhibitor, GMP Grade	038410
GMP-M036	Pyrophosphatase, Inorganic (yeast), GMP Grade	038438
GMP-E127	DNase I, GMP Grade	038364
GMP-M062	Vaccinia Capping Enzyme, GMP Grade	038330
GMP-M072	mRNA Cap 2'-O-Methyltransferase, GMP Grade	038439
GMP-E224	RNase R, GMP Grade	039011

## 5 Files Support for IND

No.	File Name	No.	File Name
1	Business License	9	Certificate of Analysis (COA)
2	IS09001 Certification	10	Manufacture Process
3	Declaration of GMP Quality Management System	11	Testing Methods (including mycoplasma, sterility, endotoxin, activity testing, etc.)
4	Declaration of TSE/BSE	12	Stability Report
5	Material Safety Data Sheet (MSDS)	13	Key Materials List
6	Certificate of Origin (COO)	14	Verification Report of Cell Bank
7	Equipments List of Production and Inspection	15	Verification Report on Virus Removal Process
8	Quality Standards	16	Other Required Documents (negotiable)

**Distributed by:**

**CliniSciences Group**

## 6 GMP Grade Products

### (1) Recombinant anti-Human CD3 mAb/CD28 mAb, GMP Grade

The first signal for T cell activation *in vivo* comes from the T cell receptor (TCR) specifically recognizing the MHC molecule-antigen polypeptide complex on the surface of the antigen presenting cell (APC); the interaction between the co-stimulatory molecule CD28 on the surface of the T cell and its ligand B7 (CD80/86) (second signal) can enhance the activation and proliferation of T cells. The activation of T cells *in vitro* can be achieved by binding anti-CD3 mAb to CD3 molecules, and anti-CD28 mAb can bind CD28 as a costimulatory molecule for T cell activation.

Novoprotein provides GMP grade humanized Anti-CD3 mAb and Anti-CD28 mAb. Anti-CD3 mAb is derived from the OKT3 through humanized modification, targeting CD3ε; Anti-CD28 mAb can bind CD28 as a costimulatory molecule for T cell activation.




Figure 2. Schematic diagram of the mechanism of CD3/CD28 antibodies in T cell activation

### Product Information

Cat. No.	Product Name
GMP-A018	Recombinant anti-Human CD3 mAb
GMP-A063	Recombinant anti-Human CD28 mAb

Distributed by:

CliniSciences Group

## (2) NovoNectin<sup>®</sup>, GMP Grade

NovoNectin<sup>®</sup> (recombinant human Fibronectin, FN-CH296) can be used for cell attachment, spreading, differentiation, and proliferation. It can greatly improve the infection efficiency of retroviruses on mammalian cells. VLA-4 and VLA-5 on the cell surface bind to the CS-1 site and the cell-binding domain on NovoNectin<sup>®</sup>, respectively, and the retroviral vector binds to the heparin-binding domain, thus promoting the transfection efficiency of retroviruses and lentiviruses on cells. NovoNectin<sup>®</sup> can be mixed with humanized monoclonal antibody CD3 for the coating to enhance T cell amplification.

Recombinant Human NovoNectin<sup>®</sup> is expressed by *E. coli* and is manufactured using animal-origin-free raw materials. Corresponding drug products are manufactured with pharmaceutical-grade excipients. The manufacturing follows the “Good Manufacturing Practice for Drugs (2010 Revision)” with strict control of residual bacterial endotoxin, residual host cell protein, and residual exogenous DNA to ensure product quality, safety, and efficacy.

### Recommendations for Use:

#### NovoNectin<sup>®</sup> for plate coating:

1. Calculate the amount of NovoNectin<sup>®</sup> required on the basis of 5 µg/cm<sup>2</sup> coating area and dilute the protein mother liquor to 20–100 µg/mL with normal saline;
2. Add the protein mother liquor to the coated vessel to cover the vessel surface and allow it to stand at room temperature for 2 h or at 4 ° C overnight;

\* Add 0.5 mL/well to the 24-well plate; add 2 mL/well to the 6-well plate or 35 mm dish;

#### Viral transfection:

1. Add the virus supernatant onto the NovoNectin<sup>®</sup>-coated plate at 125–250 µL/cm<sup>2</sup>;

\* MOI = 1, virus titer > 1 × 10<sup>8</sup>.

2. Allow standing at 37 ° C for 4–6 h for complete virus adsorption;

\* In the case of a lower virus titer, Concentrate the virus supernatant before use.

3. Prepare the target cells to a suspension of 0.2–1 × 10<sup>5</sup> cells;


\* The target cells may be activated 24 h in advance.

4. Add the target cells to the previously processed dish at 0.5–2.5 × 10<sup>4</sup> cells/cm<sup>2</sup>;

5. Incubate at 37 ° C for 2–3 days;

\* In the case of a higher virus titer, the virus may be directly mixed with the cells and the mixture shall be added to the NovoNectin<sup>®</sup>-coated dish for culture at an appropriate density.

6. The infected cells can be cultured in a conventional way.



Measured by its ability to support Jurkat cell attachment and spreading when used as a substratum for cell culture.

## Product Information

Cat. No.

GMP-CH38

Product Name

NovoNectin<sup>®</sup>, GMP Grade

**Distributed by:**

**CliniSciences Group**




### (3) Vitronectin, GMP Grade


Vitronectin is a large glycoprotein found in blood and the extracellular matrix (ECM). Vitronectin is involved in a number of biological activities including cell adhesion, cell spreading and migration, cell proliferation, extracellular anchoring, fibrinolysis, hemostasis, and complement mediated immune defense. It is worth noting that in the culture of pluripotent stem cells, due to the clear chemical composition of the protein, Recombinant Human Vitronectin helps to reduce the instability, compared with human plasma-derived Vitronectin and basement membrane extracts (BMEs).

Recombinant Human Vitronectin is expressed by mammalian cells, and is produced with raw materials of pharmaceutical applicable level. The host protein residue, nucleic acid residue and common pathogens are strictly controlled, and the production and quality management procedures of the product comply with GMP regulations to ensure the traceability of the production process and all raw materials.

#### High Activity



Measured by its ability to support iPS cell attachment and spreading when used as a substratum for cell culture. The ED50 for this effect is 0.12µg/mL.



The activity of Recombinant Human Vitronectin (Cat. No.:GMP-C395) was compared to another commercially available product.

#### Suggested Product Dosage (Recommended Working Concentration: 0.5µg/ml)

Culture vessel	Approximate surface area	VTN
6-well plate	10 cm <sup>2</sup> /well	0.5µg/well
12-well plate	4 cm <sup>2</sup> /well	0.2µg/well
24-well plate	2 cm <sup>2</sup> /well	0.1µg/well
35-mm dish	10 cm <sup>2</sup>	0.5µg
60-mm dish	20 cm <sup>2</sup>	1µg
100-mm dish	60 cm <sup>2</sup>	3µg
T-25 flask	25 cm <sup>2</sup>	1.25µg
T-75 flask	75 cm <sup>2</sup>	3.75µg

#### Product Information

Cat. No.	Product Name
GMP-C395	Recombinant Human Vitronectin

**Distributed by:**  
**CliniSciences Group**

## (4) Cell/Organoid Culture Cytokines and Antibodies, GMP Grade

### Product Features

200+ cell therapy clients are using


Successfully assisted clients in applying for IND

High purity, high activity, and low endotoxin


High lot-to-lot consistency

Stable long-term supply

### High Activity




Measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells. The specific activity of Recombinant human IL-15 (Cat.No.:GMP-C016) is  $\geq 1.0 \times 10^7$  IU/mg, which is calibrated against human IL-15 WHO International Standard (NIBSC code: 95/554).




The activity of Recombinant human IL-15 (Cat. No.:GMP-C016) was compared to another commercially available product.

### High Lot-to-Lot Consistency



Three independent lots were tested for activity and plotted on the same graph to show lot-to-lot consistency of IL-7.



Three independent lots were tested for activity and plotted on the same graph to show lot-to-lot consistency of IL-4.

Distributed by:

CliniSciences Group

## Cell/Organoid Culture Cytokines and Antibodies

Cell Type	Cytokines & Functional Antibodies
T Cell	CD3 mAb, CD28 mAb, IL-2, IL-7, IL-15, Tscm Expender®
NK Cell	CD16 mAb, NKG2D mAb, IL-12, IL-15, IL-18, IL-21
DC Cell	IL-1 alpha, IL-1 beta, IL-4, GM-CSF, TNF alpha
MSC Cell	bFGF, BMP-4, EGF, IGF-1, IL-6, TGF-β1, PDGF-BB
iPSC-derived T Cell	bFGF, BMP-4, FLT3L, IL-3, IL-7, SCF, VEGF165
iPSC-derived NK Cell	BMP-4, FLT3L, IL-2, IL-3, IL-7, IL-15, SCF, VEGF165
iPSC-derived Macrophage Cell	bFGF, BMP-4, FLT3L, GM-CSF, IL-3, IL-7, IGF-1, M-CSF, SCF, VEGF165
iPSC-derived Dopaminergic Neurons/Neural Precursor Cell	BDNF, FGF8b, GDNF, SHH, TGF beta3
	Motoneuron: BFGF, BDNF, CNTF, GDNF, IGF-1, SHH
iPSC-derived other Nerual Cell	Astrocyte: Activin A, BMP-4, bFGF, CNTF, EGF, IGFI, NRG1Beta
	NPC/NSC: BFGF, NT-3
iPSC-derived Cardiomyocyte	Activin A, BMP-4, BFGF, DKK1, SCF, TGF-beta 1, transferrin, VEGF165, Wnt 3a
iPSC-derived Islet Cell	Activin A, KGF
Organoids	Activin A, BDNF, bFGF, EGF, FGF8b, FGF-9, GDNF, HGF, IGF-1, KGF, Noggin, OSM, PDGF-BB, R-Spondin 1, TGF-β1, VEGF165, Vitronectin

**Distributed by:**

**CliniSciences Group**

## (5) Recombinant Cas9 Nuclease, GMP Grade

### Product Features

Animal-free/ampicillin-free

High purity, > 95%

Pharmacopoeia specifications


High gene editing efficiency

Comprehensive QC testing

High stability


### Product Data

#### (1) Comparison of Cas9 Editing Efficiency in 293T Cells



The analytical results of the cell samples edited by different brands of Cas9 proteins showed that Novoprotein Cas9 protein is superior to its competitors in both editing efficiency and knockout score.

#### (2) T7EI Cleavage Efficiency



Different Cas9/sgRNA RNP complexes were delivered to the 293T cells by electroporation. After 48 h, the cells were collected and the genomic DNA was extracted for the T7EI testing. The test results showed that under the same conditions, the cleavage efficiency of Novoprotein Cas9 nuclease is superior to that of its competitors. Lane 1: marker; Lane 2: Novoprotein; Lane 3: brand A; Lane 4: brand B; Lane 5: negative control (target gene+ T7EI); Lane 6: negative control (target gene).

### Product Information

Cat. No.

Product Name

GMP-1701

Recombinant Cas9 Nuclease, GMP Grade

**Distributed by:**

**CliniSciences Group**

## (6) (High Salt Active) BenzoNuclease®, GMP Grade

BenzoNuclease® is a genetically engineered recombinant endonuclease derived from *Serratia marcescens*. BenzoNuclease® has been shown to exert a broad spectrum of substrate specificity to degrade both DNA and RNA into 5'-monophosphate-terminated oligonucleotides, which are three to five bases in length—whether single-stranded, double-stranded, linear, circular or supercoiled. It is also known as "omnipotent nuclease" because of its high efficiency in degrading any form of DNA and RNA.

BenzoNuclease® is an effective tool enzyme for the removal of all forms of DNA and RNA from biologicals, both in laboratory studies and in industrial-scale production. Through efficient nucleic acid removal, the effect and yield of follow-up experiments and production can be significantly improved, and the performance is better than other nucleic acid removal methods.

Novoprotein provides BenzoNuclease® and High Salt Active BenzoNuclease® (HSAB) produced under GMP conditions.

Reliable quality and stable supply to assist your product development and production.

### Product Applications

- Used to remove exogenous nucleic acid from vaccine products, reduce the risk of residual toxicity of nucleic acid and improve product safety
- Used to reduce the viscosity of feed liquid caused by nucleic acid, shorten the processing time and increase the protein yield
- Used to remove nucleic acid winding on the surface of particles (viruses, inclusion bodies, etc.), which is conducive to the release and purification of particles
- Used to prepare samples for ELISA, column chromatography, 2D electrophoresis and western blot analysis. The resolution and recovery can be improved after BenzoNuclease® treatment
- Used to prevent cell clumping
- .....



Viral Vector Vaccine



Recombinant Protein and Antibody Drugs



Gene Therapy



Cell Therapy

Distributed by:

CliniSciences Group

# BenzoNuclease<sup>®</sup>, GMP Grade

## Strict Quality Control and Production Standards Create Reliable BenzoNuclease<sup>®</sup>:

- No animal-derived and human-derived ingredients, no ampicillin antibiotic
- Protein purity  $\geq 99\%$
- Without protease activity
- Bacterial Endotoxin level  $< 0.01\text{EU/KU}$

Element Standard	Criteria
Appearance	clear, transparent solution
Visible Particles	meet the specification
pH	7.5-8.5
Activity	250-400U/ $\mu\text{l}$
Specific Activity	$\geq 1.1 \times 10^6 \text{U/mg}$
Purity	$\geq 99\%$
Protease Activity	no protease activity detectable
Bacterial Endotoxins	$< 0.01 \text{EU/KU}$
Host-cell Protein Residues	$\leq 0.005\%$
Sterility	meet the specification
Heavy Metal	$\leq 10 \text{ppm}$

## Product Information



Cat. No.	GMP-1707 <small>DMF Filled</small>	GMP-1709
Product Name	BenzoNuclease <sup>®</sup> , GMP Grade	BenzoNuclease <sup>®</sup> (Tag-free), GMP Grade
Package size	100KU/200KU/5000KU	100KU/500KU/5000KU
Molecular Weight	32kDa	$30 \pm 3.0\text{kDa}$
pI	6.99	6.2
Tag	6 $\times$ His	Tag-free
Purity	$\geq 99\%$ (SDS-PAGE, SEC-HPLC)	
Optimal pH	8	
Optimal Temperature	37°C	
Cofactor	1-10 mM $\text{Mg}^{2+}$	
Formulation	20 mM Tris-HCl pH8.0, 2 mM $\text{MgCl}_2$ , 20 mM NaCl, 50% (v/v) glycerol	
Storage	Store at $-20^\circ\text{C} \pm 5^\circ\text{C}$ . Avoid repeated freeze-thaw cycles. Properly stored BenzoNuclease <sup>®</sup> is stable for at least 24 months.	
Unit Definition	In a 2.625 mL reaction system at 37°C and pH 8.0, one unit of BenzoNuclease <sup>®</sup> is defined as the amount of enzyme that causes a change in absorbance at 260 nm of 1.0 absorption units within 30 minutes.	

Distributed by:

**CliniSciences Group**

# BenzoNuclease® ELISA Kit

BenzoNuclease® ELISA Kit can detect and quantitatively analyze the residue of BenzoNuclease® in viral vectors and vaccine production with high sensitivity and specificity. The sensitivity is 0.014ng/mL, and the detection range is 0.014ng/mL-10ng/mL.



- Step 1: Test sample added to the well
- Step 2: Detection Antibody added to the well
- Step 3: Colorimetric detection with TMB substrate

## Product Information

Cat. No.	Product Name	Size
PA018	BenzoNuclease® ELISA Kit	96T


# High Salt Active BenzoNuclease®, GMP Grade

High Salt Active BenzoNuclease® is a GMP Grade nuclease developed as the most efficient solution for removal of both single and double stranded DNA and RNA at high salt conditions. This nonspecific endonuclease has peak activity at 500mM salt concentrations.

For Adenoviruses and Adeno-Associated Viruses (AAVs), which are often harvested from crude cell lysate, salt is typically added to such lysates to reduce viral aggregation, facilitating more effective nuclease action to digest residual DNA. The high salt tolerance of HSAB is particularly beneficial.


## Optimal Temperature

The best temperature for HSAB to degrade nucleic acid is 37°C, and it is active in the range of 0~42°C.





## Optimal PH

The optimum reaction pH for HSAB to degrade nucleic acid is 8, and it is active in the range of 6~10.



## Effect of Monovalent Cations

HSAB has the best enzyme activity in the high salt system of 500mM, and is highly suitable for the purification needs of AAVs, Ads and other viruses.



## Product Information

Cat. No.	Product Name	Size
GMP-1711	High Salt Active BenzoNuclease® (HSAB), GMP Grade	100KU/500KU/5000KU



## (7) mRNA Raw Materials and Reagents, GMP Grade

With the continuous progress of life science, mRNA is used as a drug in the fields of disease treatment and vaccine. The mRNA synthesized in vitro in 1990 was expressed in cells for the first time, and the mRNA vaccine in 2020 played an important role in the fight against the COVID-19. mRNA-based therapeutic modalities have caused a revolution in medicine.

mRNA vaccines have the advantages of fast development speed, flexibility, simple production process, platform-based, easy to expand production capacity, can be used for precise and personalized treatment, and can present multiple antigens at one time. The structure of mRNA vaccine is simple and can be synthesized in batches, and can be screened in high throughput during the research and development stage, which significantly saves research and development time. Effective mRNA vaccines have the characteristics of high yield, strong stability, low autoimmunity, high expression, strong antigen specificity of expression, and good stimulation effect.

Novoprotein provides total solution for RNA vaccines/drugs to solve difficulties in development and production. All products are produced with pharmaceutical grade raw materials, strictly control host protein residues, nucleic acid residues, etc., in line with GMP standards of production and quality management procedures, to ensure that the production process and all raw materials can be traced.

### Advantages

#### Leader in the industry

"Two clients' vaccines" have been authorized Chinese and overseas EUA  
Promoted nearly "10 projects" into the clinical trials

#### Stable supply

Annual capacity for "5 billion doses"

#### Qualified products for export

"FDA DMF" filed  
"HALAL" certification

#### Quality assurance

"ISO" certification


#### Reliable and efficient

Complete range of products  
"Total solutions"


Data as of April 30, 2024

### Applications (Novoprotein's IVT and capping system)


	Yield(μg)	Capping efficiency	Integrity	dsRNA content(ng/μg)
eGFP	220	98.56%	94.1%	0.015
Luciferase	183	99.97%	92.1%	0.060




Capillary electrophoresis showed good integrity and high purity of the transcript.



Mass spectrometry detection shows high capping efficiency.



eGFP and luciferase mRNA were expressed successfully after transfection.



Distributed by:

CliniSciences Group

## GMP Grade mRNA Raw Materials and Reagents List

### mRNA Preparation

Application	Cat. No.	Product Name
Plasmid Linearization	GMP-RE057	BspQI, GMP Grade <small>DMF Filled</small>
	GMP-EB057	10×BspQI Reaction Buffer, GMP Grade
	GMP-RE026	BsaI, GMP Grade <small>DMF Filled</small>
	GMP-RE036	BsaI ( <i>E. coli</i> ), GMP Grade
	GMP-EB026	10×BsaI Reaction Buffer, GMP Grade
	GMP-RE015	XbaI, GMP Grade <small>DMF Filled</small>
	GMP-EB015	10×XbaI Reaction Reaction Buffer, GMP Grade
In Vtro Transcription	GMP-E121-H200	T7 RNA Polymerase, GMP Grade <small>DMF Filled</small>
	GMP-EB231	10×Transcription Buffer, GMP Grade
	GMP-E122-H200	T7 RNA Polymerase 2.0, GMP Grade
	GMP-E125	RNase Inhibitor, GMP Grade <small>DMF Filled</small>
	GMP-M036	Pyrophosphatase, Inorganic (yeast), GMP Grade <small>DMF Filled</small>
	GMP-E131	T7 RNA Transcription Enzyme Mix, GMP Grade
	GMP-S033D-S036D	NTP (200mM Tris Solution), GMP Grade
dsDNA Template Digestion	GMP-E127	DNase I, GMP Grade <small>DMF Filled</small>
mRNA Capping	GMP-M062	Vaccinia Capping Enzyme, GMP Grade <small>DMF Filled</small>
	GMP-M072	mRNA Cap 2'-O-Methyltransferase, GMP Grade <small>DMF Filled</small>
	GMP-EB62	10×Capping Reaction Buffer, GMP Grade
	GMP-S062	SAM (32mM), GMP Grade
	GMP-S024N	GTP, GMP Grade (10mM)
mRNA Tailing	GMP-M012	<i>E. coli</i> Poly(A) Polymerase, GMP Grade
	GMP-EB12	10×Poly(A) Polymerase Buffer, GMP Grade
	GMP-S023N	ATP, GMP Grade (10mM)

### circRNA Preparation

Application	Cat. No.	Product Name
circRNA Preparation and Purification	GMP-M050	T4 RNA Ligase 2, GMP Grade
	GMP-E224	RNase R, GMP Grade <small>DMF Filled</small>
	GMP-EB224	10×RNase R Buffer, GMP Grade

Distributed by:

**CliniSciences Group**

## Purification

Application	Cat. No.	Product Name
mRNA Purification	N243	RNA Clean Beads
	S125	Lithium Chloride Precipitation Solution

## mRNA Substance Quality Control

Application	Cat. No.	Product Name
mRNA Capping Detection	CD001	mRNA Capping Detection Sample Preparation Kit (Beads)
	CD002	mRNA Capping Detection FlashPrep Kit
	E124	RNase H
	E134	Thermostable RNase H
mRNA Tailing Detection	E151	RNase T1
	E242	NovoNGS® mRNA Magnetic Isolation Kit
mRNA Raw Material Enzyme Residue Detection	PA101	Pyrophosphatase, Inorganic ELISA Kit
	PA102	T7 RNA Polymerase ELISA Kit
	PA105	RNase Inhibitor ELISA Kit
dsRNA Content Detection	RD017	NovoFast dsRNA ELISA Kit
RNase Residue Detection	DT007	RNase Detection Kit
DNase Residue Detection	DT009	DNase Detection Kit
DNA Template Residue Detection	E106	NovoStart® Probe qPCR SuperMix (UDG)
	E406	NovoStart® High-Specificity Probe qPCR SuperMix (UDG)
<i>E. coli</i> HCD Detection	DR001	NovoStart® <i>E. coli</i> DNA Residue Detection Kit
mRNA Enzymes Identification	PA007	mRNA Enzymes DIBA Kit

## Catalog mRNA

Application	Cat. No.	Product Name
Reporter Gene/ Functional Gene mRNA	MR008 / MR010	eGFP mRNA / eGFP mRNA (N1-Me-Pseudo UTP)
	MR009 / MR011	Luciferase mRNA / Firefly Luciferase mRNA (N1-Me-Pseudo UTP)
	MR201	eGFP circRNA
	MR202	Luciferase circRNA
	MR105	mCherry mRNA (N1-Me-Pseudo UTP)
	MR015	OVA mRNA (N1-Me-Pseudo UTP)
	MR016	hEPO mRNA (N1-Me-Pseudo UTP)
	MR107 / MR019	Cas9 mRNA / Cas9 mRNA (N1-Me-Pseudo UTP)
	GMP-MR005	piggyBac mRNA, GMP Grade

Distributed by:

**CliniSciences Group**

# 7 GMP Grade Products Summary

## CGT GMP Grade Relative Products

Cat. No.	Product Name
GMP-B038	anti-Human CD3/CD28 Beads
GMP-A037	Recombinant anti-Human 4-1BB mAb
GMP-A091	Recombinant anti-Human CD16 mAb
GMP-A063	Recombinant anti-Human CD28 mAb
GMP-A018	Recombinant anti-Human CD3 mAb
GMP-A052	Recombinant anti-Human CD52 mAb
GMP-A065	Recombinant anti-Human CD56 mAb
GMP-A075	Recombinant anti-Human NKG2D mAb
GMP-C687	Recombinant Human Activin A
GMP-C076	Recombinant Human BDNF
GMP-C029	Recombinant Human EGF
GMP-C046	Recombinant Human FGF basic
GMP-C798	Recombinant Human FGF-8b
GMP-C198	Recombinant Human FGF-9
GMP-CA82	Recombinant Human Flt-3 Ligand
GMP-C226	Recombinant Human GDNF
GMP-CC79	Recombinant Human GM-CSF
GMP-C414	Recombinant Human GPC3
GMP-CJ72	Recombinant Human HGF
GMP-CI57	Recombinant Human IFN-gamma
GMP-C070	Recombinant Human IL-1 alpha
GMP-CG93	Recombinant Human IL-1 beta
GMP-C013	Recombinant Human IL-2 <small>DMF Filed</small>
GMP-CD66	Recombinant Human IL-2
GMP-CF63	Recombinant Human IL-3
GMP-CD03	Recombinant Human IL-4
GMP-C009	Recombinant Human IL-6
GMP-CD47	Recombinant Human IL-7 <small>DMF Filed</small>
GMP-CI58	Recombinant Human IL-12
GMP-C016	Recombinant Human IL-15 <small>DMF Filed</small>
GMP-CD72	Recombinant Human IL-18
GMP-CC45	Recombinant Human IL-21
GMP-CH73	Recombinant Human KGF
GMP-C023	Recombinant Human LR3 IGF-1

Distributed by:

**CliniSciences Group**

Cat. No.	Product Name
GMP-CG11	Recombinant Human MICA
GMP-CB89	Recombinant Human Noggin
GMP-C099	Recombinant Human OSM
GMP-C199	Recombinant Human PDGF-BB
GMP-CX83	Recombinant Human R-Spondin 1
GMP-CI56	Recombinant Human sCD40 Ligand
GMP-CD53	Recombinant Human SCF
GMP-CA59	Recombinant Human TGF-beta 1
GMP-C008	Recombinant Human TNF-alpha
GMP-CJ95	Recombinant Human TPO
GMP-CR96	Recombinant Human VEGF165
GMP-C395	Recombinant Human Vitronectin
GMP-1701	Recombinant Cas9 Nuclease, GMP Grade <small>DMF Filled</small>
GMP-1707	BenzoNuclease®, GMP Grade <small>DMF Filled</small>
GMP-1709	BenzoNuclease® (Tag-free), GMP Grade
GMP-1711	High Salt Active BenzoNuclease® (HSAB), GMP Grade
GMP-CH38	NovoNectin® GMP Grade
GMP-1647	Tscm Expander®

## | mRNA Vaccine GMP Grade Relative Products

Turn to P15-P16 for the details.

# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mausier-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraaijenhoffstraat 137A,  
1018RG Amsterdam, - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences AS  
Address: c/o MerVerdi Munkerudtunet 10  
1164 Oslo - Norway  
Telephone: +47 21 988 882  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 916 750 700  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences AG  
Address: Fracht Ost Flughafen Kloten  
CH-8058 Zürich - Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
IWeb: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

