


General information

The **LEctPROFILE®** kits allow **efficient** evaluation of crude or purified glycoconjugate interactions (*i.e.* synthetic molecules or glycoconjugates) with lectins by a simple measurement of absorbance or fluorescence. The LEctPROFILE® kits enable a **fast measurement** (below 3 h) and are **easily accessible to all**.

Each kit is composed of a 96-well plate containing immobilised target lectin, and supplemented with the corresponding tracer and the revealing solution.

All our lectins are control quality proofed by a standardized procedure assessed by interactions with specific neoglycoproteins or glycoproteins through GLYcoPROFILE® method.



LEctPROFILE® kits are designed for relevant applications such as:

- 1) Detection of desirable or undesirable glycans during the early stage of quality control of biotherapeutics;
- 2) Identification of glyco-biomarkers related to cells behaviour and/or pathology;
- 3) Analysis and characterisation of the glycosylation level of glycoproteins contained in biological fluids or to target the sialylation status of relevant serum biomarkers;
- 4) Research of antagonist for lectins of interest that are well-known to be involved in particular disorder or in biological mechanism.


Distributed by:




How to use LEctPROFILE® kit?

The evaluation of compound interactions with lectins on LEctPROFILE® kits are possible by two different modes:

- **Direct Binding** is used to evaluate the potential interaction of compounds ranging from pure molecule to complex mixtures (glycocojugate(s), complex carbohydrates or glycomimetics). Previous labeling of target molecule(s) by biotinylation or by fluoresceinylation is required for readout.



- **Competition Mode**: Competitive interaction between a sample and a specific labelled tracer known to have a good affinity for the lectin without preliminary labeling of sample. This mode is typically used for the evaluation of the IC50-value (*i.e.* concentration corresponding to 50% of lectin inhibition), for the screening of potential candidate, for avidity comparison or for batch to batch monitoring.



Distributed by:



List of LEctPROFILE® kit

Reference	Name	Specificity	Kit Content	Analysis mode
LK01	FimH LEctPROFILE® kit	Mannosylated glycans structures	- 1 x 96 microplate well for fluorescence detection - FimH tracer solution - Streptavidine-DTAF solution	Direct Binding: Up to 10 samples analysed in triplicate at 3 concentrations Competition Mode: Up to 6 to 8 samples analysed in duplicate
LK02	DC-SIGN & Langerin LEctPROFILE® kit	Mannosylated glycans structures	- 1 x 96 microplate well for fluorescence detection - DC-SIGN&Langerin tracer solution - Streptavidine-DTAF solution	Direct Binding: Up to 7 samples analysed in triplicate at 2 concentrations Competition Mode: Up to 3 samples analysed in duplicate
LK03	Galectin 3 LEctPROFILE® kit	<i>B</i> -Galactoside glycans structures	 1 x 96 microplate well for absorbance detection Galectin-3 tracer solution Streptavidine-Peroxydase solution OPD Tablets 	Direct Binding: Up to 10 samples analysed in triplicate at 3 concentrations Competition Mode: Up to 6 to 8 samples analysed in duplicate
LK04E LK04G LK04M	Cells LEctPROFILE® kit (E: eucaryotic cells, G: germinal cells and M: Microbial cells)	Man, GalNAc, Galβ/α, GlcNAc, Fuc, Sialic acid, Rha	- 3 x 96 microplate well for fluorescence detection - CFSE Labelling solution	3 differents cells types analysed in triplicate on 22 lectins
LK05	α2,3_α2,6 LEctPROFILE® kit	Sialic acid (α2,3/α2,6) glycans structures	- 1 x 96 microplate well for fluorescence detection - α2,3 tracer solution - α2,6 tracer solution - Streptavidine-DTAF solution	Direct Binding: Up to 7 samples analysed in triplicate at 2 concentrations Competition Mode: Up to 3 samples analysed in duplicate
LK06	Galili LEctPROFILE® kit	Galα1,3Gal glycans structures	1 x 96 microplate well for absorbance detection Galili tracer solution Streptavidine-Peroxydase solution OPD Tablets	Direct Binding: Up to 10 samples analysed in triplicate at 3 concentrations Competition Mode: Up to 6 to 8 samples analysed in duplicate
LK07	ASGPR LEctPROFILE® kit	GalNAc glycans structures	- 1 x 96 microplate well for fluorescence detection - ASGPR tracer solution - Streptavidine-DTAF solution	Direct Binding: Up to 10 samples analysed in triplicate at 3 concentrations Competition Mode: Up to 6 to 8 samples analysed in duplicate
LK08	Glycoprotein LEctPROFILE® kit	Man, GalNAc, Galβ/α, GlcNAc, Fuc, Sialic acid, Rha, Xyl	- 3 x 96 microplate well for fluorescence detection - Biotin Labelling solution - Streptavidine-Peroxydase solution - Dialysis units	3 glycoproteins at 1 concentration or 1 glycoprotein at 3 concentrations in triplicate on 24 lectins

Table 1. Description of LEctPROFILE® kits.

Distributed by:





Additionally, and under request, the customization and development of a LEctPROFILE® kit according to your needs is possible. Indeed, we have access to a large choice of lectin from natural or recombinant sources (bacterial and human) that have been already validated in our GLYcoPROFILE® plateform.

Stability

Reagents included in the kit are stable for 6 months at +4°C. LEctPROFILE plate is stable for 6 months at -20°C.

Intended use: The LEctPROFILE kits are for research use only

Distributed by:

CliniSciences Group

Austria

Company: CliniSciences GmbH Address: Sternwartestrasse 76, A-1180

Wien - Austria

Telephone: +43 720 115 580 Fax: +43 720 115 577

Email: oesterreich@clinisciences.com
Web: https://www.clinisciences.com



Belgium

Company: CliniSciences S.R.L Address: Avenue Stalingrad 52, 1000

Brussels - Belgium Telephone: +32 2 31 50 800 Fax: +32 2 31 50 801

Email: <u>belgium@clinisciences.com</u> Web: <u>https://www.clinisciences.com</u>



Denmark

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: danmark@clinisciences.com Web: https://www.clinisciences.com



Finland

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fay: +45 89 884 064

Email: suomi@clinisciences.com
Web: https://www.clinisciences.com



France

Company: CliniSciences S.A.S Address: 74 Rue des Suisses, 92000

Nanterre- France Telephone: +33 9 77 40 09 09

Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: https://www.clinisciences.com



German

Company: Biotrend Chemikalien GmbH Address: Wilhelm-Mauser-Str. 41-43, 50827 Köln - Germany

Telephone: +49 221 9498 320 Fax: +49 221 9498 325 Email: info@biotrend.com Web: https://www.biotrend.com



Icelan

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: <u>island@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>



Ireland

Company: CliniSciences Limited Address: Ground Floor, 71 lower Baggot street Dublin D02 P593 - Ireland Telephone: +353 1 6971 146

Fax: +353 1 6971 147

Company: CliniSciences AS

Telephone: +47 21 988 882

Email: norge@clinisciences.com

Web: https://www.clinisciences.com

1164 Oslo - Norway

Email: <u>ireland@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>

Address: c/o MerVerdi Munkerudtunet 10



Italy

Company: CliniSciences S.r.I Address: Via Maremmana inferiore 378 Roma 00012 Guidonia Montecelio - Italy Telephone: +39 06 94 80 56 71 Fax: +39 06 94 80 00 21

Email: italia@clinisciences.com
Web: https://www.clinisciences.com



lothorlando

Company: CliniSciences B.V. Address: Kraijenhoffstraat 137A, 1018RG Amsterdam, - Netherlands Telephone: +31 85 2082 351 Fax: +31 85 2082 353

Email: nederland@clinisciences.com Web: https://www.clinisciences.com



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532

Email: polska@clinisciences.com
Web: https://www.clinisciences.com



Portugal

Company: Quimigen Unipessoal LDA Address: Rua Almada Negreiros, Lote 5, Loja 14, 2615-275 Alverca Do Ribatejo - Portugal Telephone: +351 30 8808 050

Fax: +351 30 8808 052 Email: info@quimigen.com Web: https://www.quimigen.pt



Spair

Company: CliniSciences Lab Solutions Address: C/ Hermanos del Moral 13 (Bajo E), 28019, Madrid - Spain Telephone: +34 916 750 700 Fax: +34 91 269 40 74

Email: <u>espana@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>



Sweden

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: sverige@clinisciences.com
Web: https://www.clinisciences.com



Switzerland

Company: CliniSciences AG Address: Fracht Ost Flughafen Kloten CH-8058 Zürich - Switzerland Telephone: +41 (044) 805 76 81 Fax: +41 (044) 805 76 75

Email: <u>switzerland@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511

or +44 (0) 330 684 0982 Fax: +44 (0)1753 208 899 Email: <u>uk@clinisciences.com</u> IWeb: <u>https://www.clinisciences.com</u>



USA

Company: Biotrend Chemicals LLC Address: c/o Carr Riggs Ingram, 500 Grand Boulevard, Suite 210 Miramar Beach, FL 32550- USA

Telephone: +1 850 650 7790 Fax: +1 850 650 4383 Email: <u>info@biotrend-usa.com</u>

Web: https://www.biotrend-usa.com

