

Distributed by:

CliniSciences Group

ATLAS ANTIBODIES

MolBoolean™

BEYOND PROXIMITY LIGATION

UNLOCK THE SECRETS OF PROTEIN-PROTEIN INTERACTIONS IN CELLS AND TISSUE

PRODUCT DESCRIPTION

MolBoolean™ mouse/rabbit is a kit for in situ protein proximity analysis in tissue and cells.

The MolBoolean™ assay utilizes a proprietary oligonucleotide setup that enables the simultaneous detection of both free and interacting (~40 nm proximity) fractions for two protein targets (protein A, protein B and interaction proteins AB).

MolBoolean™ Mouse/Rabbit should be performed with the user's primary antibodies of choice (raised in mouse and rabbit). MolBoolean™ relies on anti-mouse and anti-rabbit secondary proximity probes and rolling circle amplification (RCA) as a mean to amplify signal.

A series of molecular steps, performed in the MolBoolean™ assay, incorporates information into the amplified products indicating antibody target engagement of individual versus interacting proteins. This information is converted to fluorescent signals by the binding of detection reporters, allowing visualization on a conventional fluorescent microscope.

An image analysis software is provided to segment and differentiate fluorescent signals, enabling the relative quantification of free versus interacting fractions for the two analyzed protein targets in the sample.



MolBoolean™ is based on the 'boolean logic' for the simultaneous detection of free and interacting protein fractions.

WHY MOLBOOLEAN?

- Complete spatial quantitative analysis of protein-protein interactions by **simultaneous detection** of free and interacting proteins.
- Accurate **quantification by normalization** of interaction data to total target protein levels.
- Biologically relevant data **without the need** for engineered protein expression.
- **1000-fold** increased fluorescence signal by Rolling Circle amplification, allowing detection and quantification of low abundant proteins.
- **Adaptable** to different research needs.
- **Universal** kit that can be used with the customer's choice of primary antibodies.
- **Validated** in both cells and tissue.



The MolBoolean™ kit contains 15 separate tubes. The kit volume is 4.8 ml and covers approximately 120 assays in cells (40 µl/assay) and 60 assays in tissue (80 µl/assay).

Distributed by:

CliniSciences Group

MOLBOOLEAN: OVERCOMING IN SITU PLA LIMITATIONS


MolBoolean™ offers several advantages over traditional methods like in situ PLA, addressing key gaps in data interpretation and providing more reliable results:

- **Normalization of Data:** With MolBoolean™, the number of protein interactions can be normalized to the total number of target proteins. This is crucial because target protein levels can be influenced by various factors such as cell treatments or disease states. By providing a normalized measure, MolBoolean™ allows for more accurate comparisons between samples.
- **Detection of Interacting and Non-interacting Fractions:** Unlike traditional methods, MolBoolean™ can detect both interacting and non-interacting endogenous fractions for two protein targets in cells and tissues. The more comprehensive detection capability increases the spatial information that can be extracted from the experiment, allowing researcher to determine locations of positive or negative protein interactions.
- **Consistent Molecular Process:** The Rolling circle products (RCPs) stemming from free or interacting proteins are generated through the same molecular process steps. This ensures uniform signal efficiencies, reducing uncertainty during data analysis.

Overall, MolBoolean™ offers a more reliable and comprehensive solution for detecting protein interactions, addressing key gaps in data interpretation, and providing researchers with more accurate insights into molecular interactions in cells and tissues.

Assay	Individual Proteins	Protein Interaction
IHC/ICC	✓	✗
in situ PLA	✗	✓
MolBoolean	✓	✓

MolBoolean™ compared to common protein detection assays.



E-cadherin/ β-catenin MolBoolean staining in MCF7 cells.

MolBoolean™ analysis of the interaction between E-cadherin (magenta) and β-catenin (green) in MCF7 cells, using the monoclonal anti-CDH1 (Cat. AMAb90862) and the polyclonal anti-CTNNB1 (Cat. HPA029159) antibodies from Atlas Antibodies AB. Image shows the relative quantification of free versus interacting protein fractions, indicated by the detection of rolling circle products (RCPs) in either one or two fluorescent channels: 32% free E-cadherin (magenta), 19% free β-catenin (green) and 49% E-cadherin/β-catenin complex (white). Data is normalized to total target protein levels (total RCPs).

ATLAS ANTIBODIES

MolBoolean™ is a trademark of Atlas Antibodies AB

Product is for research use only. ©Atlas Antibodies AB 2024.

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>