



PRIME-XV MSC XSFM Dual Component Kit with or without Phenol Red

SIMPLIFIED AND EFFICIENT LARGE-SCALE MSC EXPANSION

Distributed by:

CliniSciences Group

Maximize cell expansion and differentiation within a fully closed system

Mesenchymal stromal/stem cells' (MSCs) proliferation, self-renewal, and multipotent differentiation potential make them cornerstones of regenerative medicine and cell therapy development. Because MSCs constitute a small percentage of their source tissues and isolation processes show variable efficiency, the ability to expand these progenitor cells *in vitro* can limit therapeutic efficacy.

Deliver greater numbers of MSCs by moving to an automated closed-system platform

- Maintains lot-to-lot consistency with xeno- and serum-free media
- Enables up to an 8-fold MSC expansion (donor-dependent)
- Yields 200×10^6 cells in 7 days of expansion with viability at or above 90%
- Does not contain antibiotics, with or without phenol red
- Extensive QC testing on all lots including:
 - Mycoplasma (USP <63>), functionality, sterility (USP <71>), and endotoxin (USP <85>)



PRIME-XV MSC XFSM Dual Component Kit with or without Phenol Red

Our closed-system solutions include basal media and supplement in a complete kit for use with Quantum or other compatible hollow fiber bioreactors. A simplified and efficient GMP solution for large-scale MSC expansion.

- Simplifies and expedites the process of MSC expansion for large-scale solutions
- Minimizes risk of contamination and cell stress, and increases efficiency
- Provides a fully closed system solution for Quantum* or other compatible hollow fiber bioreactors

PRIME-XV MSC XFSM Dual Component Kit is comprised of two components; PRIME-XV MSC XFSM Base media and PRIME-XV MSC XFSM Supplement.

8-fold
expansion
in xeno- and
serum-free
conditions

*Quantum Bioreactor protocol available in IFU (Product Insert)

Distributed by:
CliniSciences Group

Consistent, Reproducible MSC Expansion

PRIME-XV MSC XSFM Dual Component Kit with or without Phenol Red is formulated to deliver expansion and viability of MSCs over several passages in culture, while maintaining multipotency and functional characteristics.

Results from three different donors demonstrate some donor-to-donor variability. This supports a yield of 200×10^6 total MSCs from an initial seed of 25×10^6 cells in 7 days of expansion and helps maintain stemness of ADMSCs post expansion.

Large-scale expansion of MSCs in a closed system format




Figure 1. The dual bag format keeps the critical, temperature-sensitive supplement frozen until ready to use. The supplement and basal bag are aseptically welded and mixed to create the complete medium.

- A)** 10 mL PRIME-XV MSC XSFM Supplement stored at -20°C
- B)** 5 L PRIME-XV MSC XSFM Base with Phenol Red or PRIME-XV MSC XSFM Base without Phenol Red in bag stored at $2 - 8^{\circ}\text{C}$
- C)** Quantum Bioreactor or another compatible bioreactor system




Figure 2. Results from the 7-day expansion PRIME-XV MSC XSFM Dual Component Kit with Phenol Red show that viability is consistently maintained at or above 90%.




Figure 3. ADMSCs maintain their stem cell phenotype at Day 7 post-expansion in the Terumo Quantum bioreactor with PRIME-XV MSC XSFM Dual Component with Phenol Red and in PRIME-XV MSC Expansion XSFM (Catalog #91149). The expression of negative markers CD34 and CD45 remained well below 2% over the course of the culture.

Distributed by:

CliniSciences Group

Ordering Information

Product Description	Catalog #	Size*	Additional Information
PRIME-XV MSC XSFM Dual Component Kit with Phenol Red	91149DC	10 mL Supplement 5 L Basal	Xeno-free, serum-free, medium for MSC expansion.
PRIME-XV MSC XSFM Dual Component Kit without Phenol Red	91214DC	10 mL Supplement 5 L Basal	Xeno-free, serum-free, medium for MSC expansion.

Related Products

Product Description	Catalog #	Size*	Additional Information
PRIME-XV MSC Expansion XSFM	91149	250 mL 1 L	Xeno-free, serum-free, medium for MSC expansion.
PRIME-XV FreezIS DMSO-Free	91140	10 mL 100 mL	Protein-free, chemically defined, animal component-free cryopreservation medium. Does not contain DMSO.
PRIME-XV FreezIS	91139	10 mL 100 mL	Protein-free, chemically defined, animal component-free cryopreservation medium. Contains DMSO.
PRIME-XV Adipogenic Differentiation SFM	91137	100 mL	Serum-free adipogenic differentiation medium.
PRIME-XV Chondrogenic Differentiation XSFM	91138	100 mL	Xeno-free, serum-free chondrogenic differentiation medium.
PRIME-XV Human Fibronectin	31002	1 mg	Human plasma-derived fibronectin, carrier-free.

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oessterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: danmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkeordtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

