
# κλιτυς


# VLP & Nanodisc Displayed Antigens

Multi-pass transmembrane proteins displayed on virus-like particles or copolymer nanodiscs



### VLP-Displayed Multi-pass Transmembrane Proteins

Virus-like particles (VLPs) originate from the capsid proteins of viruses. They are tiny nanoparticles formed by the self-assembly of one or more types of capsid proteins. VLPs do not contain the infectious viral genome, making them relatively safe during production operations. With technological advancements, VLPs have demonstrated significant advantages in presenting antigens for the discovery of antibodies against challenging targets due to their effectiveness in efficiently activating both humoral and cellular immune responses in the body.



We have successfully developed a series of full-length multi-pass transmembrane proteins displayed on VLPs. These targets include Claudin 18.2, GPRC5D, Claudin 6, CD20, and more. Our VLP-displayed proteins overcome the technical barriers of membrane protein expression and stability. Additionally, the VLP structure increases immunogenicity of the target antigen for immunization campaigns. Our VLP proteins have undergone rigorous biological activity verification including SPR, BLI, and ELISA. These products feature high sensitivity, high specificity, and high stability.

### VLP | Product Advantages



Full-length sequence expression with native structure and conformation



Mammalian expression



Boosted immunogenicity



Activity Verification via ELISA, SPR, FACS, etc.

### VLP | Product Applications





Animal Immunization In vitro efficacy and Pharmacokinetic Studies





Antibody Screening

ELISA, SPR, BLI and other analytical tests

### VLP | Product Performance Validation





Figure 2. ELISA data confirms that Human Claudin 18.2 VLP can bind to Anti-Claudin 18.2 antibody with an EC50 value of 9.8ng/mL.

Catalog No. CLD-HE1822



BLI | Biotinylated Claudin 18.2 VLP

Figure 4. BLI data demonstrates high binding affinity (KD=5.73E-10 M) of Biotinylated Human Claudin 18.2 VLPs to streptavidin-labeled probes.

Catalog No. CLD-HE1822B

SPR | Biotinylated Claudin 18.2 VLP RU 25 20 15 Response 10 5 0 -5 -100 0 100 200 300 400 500 Time s

Figure 3. SPR data confirms that Biotinylated Human Claudin 18.2 VLP can specifically bind to Anti-Claudin 18.2 antibody with an affinity constant of 1.28 nM.

Catalog No. CLD-HE1822B





Figure 5. SPR data confirms that Biotinylated Human GPRC5D VLP can specifically bind to Anti-GPRC5D antibody with an affinity constant of 0.30 nM.

Catalog No. GPR-HM05PB

### Nanodisc-Displayed Multi-pass Transmembrane Proteins

KACTUS Copolymer Nanodiscs are a novel technology displaying membrane proteins in native conformations. Traditional membrane scaffold protein nanodiscs rely on detergents for membrane protein extraction and stabilization for reconstitution. Our innovative approach allows copolymer nanodiscs to be prepared directly from cells with detergent-free extraction. Our Nanodisc products offer mammalian-expressed full-length solubilized membrane proteins for analytical assays including ELISA, SPR, and BLI.



### Nanodisc | Product Advantages



Native membrane protein structure and conformation



Water-solubilized full-length proteins



Detergent-free extraction to preserve protein activity



Activity verification via ELISA, SPR, FACS

### Nanodisc | Applications



### Nanodisc | Product Performance Validation




Figure 8. ELISA data demonstrates Human GPRC5D Nanodisc can bind to Anti-GPRC5D antibody with an EC50 of 4.9ng/mL.

Catalog No. GPR-HM15P




Figure 9. Human CD3E&CD3D was immobilized on the plate, followed by Anti-Human CD3×GPRC5D Bispecific Antibody binding, followed by Biotinylated Human GPRC5D Nanodisc. Results demonstrate quality performance in ELISA to of the nanodisc to Anti-Human CD3×GPRC5D bispecific antibody with an EC50 of 0.28 µg/mL.

Catalog No. GPR-HM45PB




Figure 10. BLI data demonstrates Biotinylated Human GPRC5D Nanodisc can bind to Anti-GPRC5D Antibody with high affinity (KD=1.16nM).

Catalog No. GPR-HM45PB




Figure 11. SPR data demonstrates Human GPRC5D Nanodisc can bind to Anti-GPRC5D Antibody with a high affinity (KD = 1.47nm).

Catalog No. GPR-HM15P

# Distributed by: CliniSciences Group

### BLI | Biotinylated GPRC5D Nanodisc

ELISA | Biotinylated GPRC5D Nanodisc

### Application of GPRC5D Copolymer Nanodisc in Isolating GPRC5D Binders for Yeast Display Antibody Discovery

Our Copolymer Nanodisc effectively facilitated the identification of high-binding anti-GPRC5D antibody clones via yeast display. The Biotinylated GPRC5D Nanodisc (Cat No. GPR-HM45PB) showed strong and specific binding interactions with various monoclonal antibodies (mAb) displayed on yeast cells, as evidenced by significant fluorescence signals across multiple antibody clones at different GPRC5D membrane protein concentrations. Importantly, the GPRC5D Nanodisc displayed no binding to control yeast cells, underscoring its suitability for selective yeast display screening.




Figure 12. Plate view for Biotinylated GPRC5D Nanodisc binding various mAb displayed on yeast. Most mAb clones showed significant binding to the Biotinylated GPRC5D Nanodiscs and negative controls showed minimal or no binding, indicating specificity of the antibody clones for GPRC5D.




Figure 13. Biotinylated GPRC5D Nanodisc binding Rabbit Anti-GPRC5D Clone 1 Fab. Results show Rabbit Anti-GPRC5D Clone 1 Fab displayed on yeast cells binds with high specificity to Biotinylated GPRC5D Nanodisc.

We offer innovative customized solutions for studying difficult-to-express target antigens displayed on VLPs or Nanodiscs. These platforms provide a stable and controlled environment, enabling detailed analysis of protein structure, function, and interactions. Our custom services include full-length membrane protein expression on your choice of display platform (VLP or Nanodisc), choice of biotinylation or fluorescent labeling, choice of expression system, bulk quantities, and optional SPR analytical services.

### Service Features:



8 weeks turnaround time




HEK293 Expression



Full-Length Protein

# CXCR4 VLP Immunization Dosage & Antiserum Generation

Our team has expressed full-length CXCR4, a complex GPCR protein characterized by its seven transmembrane domains, on VLPs and demonstrated successful antiserum generation even with low dosages of CXCR4 VLP (0.5µg).



**Figure 14.** An immunization dose-range study showing as low as 0.5µg VLP can produce a strong immune response.

### **Areas of Expertise**

GPCRs			Solute Carriers		4HB/Tetraspanins	
CCR2b	CB1	GPRC5D	SLC1A4	SLC16A3	Claudin 1	Claudin 18.2
CCR4	CB2	A2AR	SLC1A5	SLC34A2	Claudin 3	CD20
CCR5	CXCR4	BILF1	SLC7A1	SLC40A1	Claudin 6	TM4SF1
CCR6	GPR75	MRGPRX2	SLC15A4	SLC44A4	Claudin 9	
CCR7	GPR77	SSTR2	SLC13A5	SLC59A1		
CCR8	GLP-1R	GCGR	SLC16A1			

### **Product List**

Catalog No.	Product	Species	Tag	Express System
A2R-HM1N1	A2AR Nanodisc	Human	C-His	HEK293
CNR-HM001	Cannabinoid receptor 1 VLP	Human		HEK293
CCR-HM02BB	Biotinylated CCR2b VLP	Human		HEK293
CCR-HM02B	CCR2b VLP	Human		HEK293
CR4-HM1N122	FITC-equivalent CCR4 Nanodisc	Human	C-His	HEK293
CCR-HM107	CCR7 Nanodisc	Human	C-His	HEK293
CD33-HM1N144	CD133 Nanodisc	Human	C-His	HEK293
CD2-HM122	CD20 VLP	Human		HEK293
CD2-HM1N37	CD20 Nanodisc Nanodisc	Human	C-His	HEK293
CD2-CM124V	CD24 VLP	Cynomolgus		HEK293
CD2-HM124V	CD24 VLP	Human		HEK293
CLD-HM40NB	Biotinylated Claudin 18.2 Nanodisc	Human	N-His,C-Avi	HEK293
CLD-HE1822B	Biotinylated Claudin 18.2 VLP	Human		HEK293
CLD-HM182FB	FITC-equivalent Biotinylated Claudin 18.2 VLP	Human		HEK293
CLD-HM10N	FITC-equivalent Claudin 18.2 Nanodisc	Human	N-His	HEK293
CLD-HM0P9	FITC-equivalent Claudin 18.2 VLP	Human		HEK293
CLD-HE1822	Claudin 18.2 VLP	Human		HEK293
CLD-HM006B	Biotinylated Claudin 6 VLP	Human		HEK293
CLD-CM006	Claudin 6 VLP	Cynomolgus		HEK293
CLD-HM006	Claudin 6 VLP	Human		HEK293
CLD-MM006	Claudin 6 VLP	Mouse		HEK293
CLD-HM009	Claudin 9 VLP	Human		HEK293
GPC-HM003	GPC3 (438-554) VLP	Human		HEK293
GPC-HE005	GPC3 VLP	Human		E.coli
GPR-HM45PB	Biotinylated GPRC5D Nanodisc	Human	C-His-Avi	HEK293
GPR-HM05PB	Biotinylated GPRC5D VLP	Human		HEK293
GPR-HM15P	GPRC5D Nanodisc	Human	C-His	HEK293
GPR-HM05P	GPRC5D VLP	Human		HEK293
GPR-CM05P	GPRC5D VLP	Cynomolgus		HEK293
GPR-MM05P	GPRC5D VLP	Mouse		HEK293
LGR-HM10N	LGR-4 Nanodisc	Human	C-His	HEK293
MR2-HM1N118	MRGPRX2 Nanodisc	Human	C-His	HEK293
SL7-HM0P29	FITC-equivalent SLC7A1 VLP	Human		HEK293
STR-HM1N1	SSTR2 Nanodisc	Human	C-His	HEK293
STR-HM002	SSTR2 VLP	Human		HEK293
TSF-HM00N	TM4SF1 Nanodisc	Human	N-His	HEK293
TSF-HM002	TM4SF1VLP	Human		HEK293
GPR-HM05CB	Biotinylated VLP Control	N/A		HEK293
CON-HM0P34	FITC-equivalent VLP Control	N/A		HEK293
VLP-HM00C	VLP Control	N/A		HEK293

# **CliniSciences** Group

#### Austria

Company: CliniSciences GmbH Address: Sternwartestrasse 76, A-1180 Wien - Austria Telephone: +43 720 115 580 Fax: +43 720 115 577 Email: oesterreich@clinisciences.com Web: https://www.clinisciences.com

#### Finland

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064 Email: suomi@clinisciences.com Web: https://www.clinisciences.com

#### Iceland

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064 Email: island@clinisciences.com Web: https://www.clinisciences.com

#### Netherlands

Company: CliniSciences B.V. Address: Kraijenhoffstraat 137A, 1018RG Amsterdam, Netherlands Telephone: +31 85 2082 351 Fax: +31 85 2082 353 Email: nederland@clinisciences.com Web: https://www.clinisciences.com

#### Portugal

Company: Quimigen Unipessoal LDA Address: Rua Almada Negreiros, Lote 5, Loja 14, 2615-275 Alverca Do Ribatejo - Portugal Telephone: +351 30 8808 050 Fax: +351 30 8808 052 Email: info@guimigen.com Web: https://www.quimigen.pt

#### Switzerland

Company: CliniSciences AG Address: Fracht Ost Flughafen Kloten CH-8058 Zürich - Switzerland Telephone: +41 (044) 805 76 81 Fax: +41 (044) 805 76 75 Email: switzerland@clinisciences.com Web: https://www.clinisciences.com

#### Belaium

Company: CliniSciences S.R.L Address: Avenue Stalingrad 52, 1000 Brussels - Belgium Telephone: +32 2 31 50 800 Fax: +32 2 31 50 801 Email: belgium@clinisciences.com Web: https://www.clinisciences.com

### France

Company: CliniSciences S.A.S Address: 74 Rue des Suisses, 92000 Nanterre- France Telephone: +33 9 77 40 09 09 Fax: +33 9 77 40 10 11 Email: info@clinisciences.com Web: https://www.clinisciences.com

### Ireland

Company: CliniSciences Limited Address: Ground Floor, 71 lower Baggot street Dublin D02 P593 - Ireland Telephone: +353 1 6971 146 Fax: +353 1 6971 147 Email: ireland@clinisciences.com Web: https://www.clinisciences.com

#### Norway

Company: CliniSciences AS Address: c/o MerVerdi Munkerudtunet 10 1164 Oslo - Norway Telephone: +47 21 988 882 Email: norge@clinisciences.com Web: https://www.clinisciences.com

Spain

Company: CliniSciences Lab Solutions Address: C/ Hermanos del Moral 13 (Bajo E), 28019, Madrid - Spain Telephone: +34 91 269 40 65 Fax: +34 91 269 40 74 Email: espana@clinisciences.com Web: https://www.clinisciences.com

### UK

Company: CliniSciences Limited Address: 11 Progress Business center, Whittle Parkway, SL1 6DQ Slough- United Kingdom Telephone: +44 (0)1753 866 511 or +44 (0) 330 684 0982 Fax: +44 (0)1753 208 899 Email: uk@clinisciences.com IWeb: https://www.clinisciences.com



#### Denmark

Germany

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064 Email: danmark@clinisciences.com Web: https://www.clinisciences.com



Company: Biotrend Chemikalien GmbH Address: Wilhelm-Mauser-Str. 41-43, 50827 Köln - Germany Telephone: +49 221 9498 320 Fax: +49 221 9498 325 Email: info@biotrend.com Web: https://www.biotrend.com

### Italy



Company: CliniSciences S.r.I Address: Via Maremmana inferiore 378 Roma 00012 Guidonia Montecelio - Italy Telephone: +39 06 94 80 56 71 Fax: +39 06 94 80 00 21 Email: italia@clinisciences.com Web: https://www.clinisciences.com

#### Poland



Company: CliniSciences sp.Z.o.o. Address: ul. Rotmistrza Witolda Pileckiego 67 lok. 200 - 02-781 Warszawa -Poland Telephone: +48 22 307 0535 Fax: +48 22 307 0532 Email: polska@clinisciences.com Web: https://www.clinisciences.com

### Sweden

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064 Email: sverige@clinisciences.com Web: https://www.clinisciences.com

### USA

Company: Biotrend Chemicals LLC Address: c/o Carr Riggs Ingram, 500 Grand Boulevard, Suite 210 Miramar Beach, FL 32550- USA Telephone: +1 850 650 7790 Fax: +1 850 650 4383 Email: info@biotrend-usa.com

Web: https://www.biotrend-usa.com

