## **Specification Sheet**

## Monoclonal Antibody Against Bovine Lactoferrin

Code	2ABLF3
Clone No.	Hyb- K1256
Manufacturing Process	Hybridoma obtained by fusing mouse myeloma cells (NS-1) with BALB/c mouse spleen cells immunized with purified lactoferrin derived from bovine colostrum is proliferated in the BALB/c mouse abdominal cavity and ascites is collected.
Specificity	This monoclonal antibody specifically reacts with bovine lactoferrin (c) and does not react with human lactoferrin.
Purification	Ammonium sulfate precipitation.
Ig Class	IgG <sub>1</sub>
Form	Saline (containing 0.1% NaN <sub>3</sub> )
Protein Concentration	0.5 mg/mL
Package	0.4 mL
Storage	2-8 (Below -20 for long term storage)
Application	ELISA

Reference



1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004 Tel : +81-3-3814-4081 Fax : +81-3-3814-5957 e-mail : info@tokumen.co.jp http://www.tokumen.co.jp/english/

This product is for *in vitro* research use only and is not intended for use humans or animals.

2ABLF31210-03