



AteloCell[®]

DATA SHEET

Native collagen Bovine dermis 5mg/mL

Cat No.: IAC-50

Lot No.: 211028

ORIGIN:	Bovine dermis
CONCENTRATION:	5.1mg/mL
SOLUTION:	1mmol/L HCl
pH:	3.12
STERILITY TEST:	Negative (medium: TGC-I and SCD)
CELL CULTURE TEST:	Normal (Human Fibroblast)
SPECIFIC OPTICAL ROTATION: $[\alpha]_D^{20}$	-412°
PURITY:	Total collagen >95% by SDS-PAGE
REFERENCES:	Refer to the AteloCell [®] website http://www.atelocell.com

FOR RESERCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.

3-14-3 Mejiro, Toshima-ku, Tokyo Japan TEL: +81-3-3950-6600 FAX: +81-3-3950-6602

<http://www.atelocell.com>