

## RB1 antibody

### Product Information

Catalog No.:	FNab07140
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	retinoblastoma 1
Synonyms:	Retinoblastoma-associated protein p105-Rb p110-RB1 pRb (Rb) pp110 RB1
Observed MW:	110 kDa
Uniprot ID :	P06400

### Application

Reactivity:	Human, Mouse, Rat, Zebrafish
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HT-1080 cells were subjected to SDS PAGE followed by western blot with FNab07140(Rb antibody) at dilution of 1:500