Ab Science, Inc.

Product Datasheet

Mouse Anti-feline CD3 ε Monoclonal Antibody


Product munber: ABS-070002

Clone: 5G-6-7-3

Size: 0.1mg

CD number: CD3

Product Form: Purified IgG – liquid



Preparation: Purified IgG prepared by affinity chromatography on Protein G from culture supernatant

Reconstitution Buffer: Reconstitute with sterile PBS containing 0.09% sodium azide

Specificity: This antibody was selected for its ability to detect feline CD3 ε molecule in flow cytometry and western blots.

Species Cross Reactivity: Do not react with canine

Immunogen: Synthetic protein sequence derived from feline CD3 ε

Ig Isotype ∶ Mouse IgG1

Application:

	Yes	No	Not	Suggested
			Determined	Working Dilution
Western Blotting	0			1/200 — 1/1000
Flow Cytometry	0			
Immunohistology			0	
ELISA	0			

Ab Science, Inc.

Optimal dilutions should be determined by each laboratory for each application.

Storage: Store at 4°C or at -20°C if preferred

References: No known references at this time

For research use only.