
TARDBP antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab08577 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 9H8 |
| IsoType: | IgG2a |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information


| | |
|--------------|--|
| Immunogen: | TAR DNA binding protein |
| Synonyms: | TAR DNA-binding protein 43 (TDP-43) TARDBP TDP43 |
| Observed MW: | 43 kDa |
| Uniprot ID : | Q13148 |

Application

| | |
|-----------------------|---|
| Reactivity: | Human |
| Tested Application: | ELISA, IP, FC, WB, IHC |
| Recommended dilution: | WB :1:5000-1:20000;IP :1:500-1:2000; IHC : 1:500-1:5000 |
| Image: | |



Immunohistochemistry of paraffin-embedded human brain using FNab08577(TARDBP antibody) at dilution of 1:1000



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab08577(TARDBP antibody) at dilution of 1:10000.