
ALB antibody

Product Information

Catalog No.:	FNab10149
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	BSA
Synonyms:	Albumin BSA ALB
Observed MW:	67 kDa
Uniprot ID :	P02769

Application

Reactivity:	Bovine
Tested Application:	ELISA