

## **Rabbit anti-Human Somatostatin Receptor SST2B (SSTR2B)**

**Cat# NB-49-022-50ul**

### **Description**

---

Rabbit anti-Human Somatostatin Receptor SST2B (SSTR2B)

### **Product Informations**

---

<b>Cat.No.</b>	NB-49-022-50ul
<b>Size</b>	50 µl
<b>Host:</b>	Rabbit
<b>Immunogen:</b>	Peptide coupled to KLH
<b>Presentation:</b>	Whole Antiserum, liquid
<b>Specificity:</b>	This antiserum recognizes the C-terminus of human SSTR2B receptor. SSTR2B is a putative splice variant of Somatostatin Receptor 2. Specificity of antibody binding was accomplished by blocking with is cognate peptide.
<b>Application:</b>	Immunohistochemistry
<b>Titer:</b>	Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.
<b>Storage:</b>	Aliquot and store at -20°C or colder. Avoid repeated freezing and thawing cycles.