

DYNC1H1 antibody

Product Information

Catalog No.:	FNab02582
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	dynein, cytoplasmic 1, heavy chain 1
Synonyms:	Cytoplasmic dynein 1 heavy chain 1 Cytoplasmic dynein heavy chain 1 Dynein heavy chain, cytosolic DYNC1H1 DHC1 DNCH1 DNCL DNECL DYHC KIAA0325
Observed MW:	532 kDa
Uniprot ID :	Q14204

Application

Reactivity:	Human, Rat, Zebrafish, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02582(KIAA0325 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02582(KIAA0325 antibody) at dilution of 1:500