

## NR2C2 antibody

### Product Information

Catalog No.:	FNab08492
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitogen-activated protein kinase kinase kinase 7
Synonyms:	Nuclear receptor subfamily 2 group C member 2 Orphan nuclear receptor TAK1 Orphan nuclear receptor TR4 Testicular receptor 4 NR2C2 TAK1 TR4
Observed MW:	64-67 kDa, 75-80 kDa
Uniprot ID :	P49116

### Application

Reactivity:	Human, Hamster, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using FNab08492(MAP3K7 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



CHO cells were subjected to SDS PAGE followed by western blot with FNab08492(MAP3K7 antibody) at dilution of 1:1000