

SFTPD antibody

Product Information

Catalog No.:	FNab08401
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	surfactant protein D
Synonyms:	Pulmonary surfactant-associated protein D (PSP-D, SP-D) Collectin-7 Lung surfactant protein D SFTPD COLEC7 PSPD SFTP4
Observed MW:	37 kDa
Uniprot ID :	P35247

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IP
Recommended dilution:	IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human lung cancer using FNaB08401(Surfactant protein D antibody) at dilution of 1:50



IP Result of anti-Surfactant protein D (IP:FNaB08401, 3ug; Detection:FNaB08401 1:300) with mouse lung tissue lysate 2800ug.