



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse TSPYL4 (KIAA0721) Polyclonal Antibody, Rabbit

Code # : MKA0721AF [FAF00118], (Lot# XXX)

Gene : mouse TSPY-like 4 (mTSPYL4, mKIAA0721)

Immunogen : GX1772 (GST-fusion protein, 251 amino acids)
TGKEGEAGAAMQEKKGLQKEKKVAGGGKEETRPRAPKINCMSLEAIDQELSNVNAQADRAFL
QLERKFGRMRRLHMQRSSFIIQNIPGFVVTAFRNHPQLSPMISGQDEDMRYMINLEVEELKQ
PRVGCKFKFIFQSNPYFRNEGLVKEYERRSSGRVVSLSTPIRWHRGQEPQAHIHNRNREGNTIPS
FFNWFSDHSLLEFDRIAIEIIGELWSNPLQYYLMGDGPRRGVVRVPPRQPVESPRSFQSG

Format : Affinity Purified Rabbit IgG

Constitution : PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity : 50 µg (250 µL)

Specificity : Specific to recombinant protein GX1772. This antibody detects mTSPYL4 protein. It also recognizes human TSPYL4 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>