

MOUSE ANTI IgE  
MONOCLONAL ANTIBODY

CLONE:	29D10
QUANTITY:	0.1mg
SPECIFICITY:	Human IgE
IMMUNOGEN:	Human IgE
Type/Subtype:	IgG1 $\kappa$
APPLICATIONS:	Sandwich ELISA(This antibody can be used as a detection antibody in a human IgE ELISA combination with the monoclonal capture antibody 23F4.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

*For research use only; not for use as a diagnostic.*

Produced by:  
Japan Clinical Laboratories, Inc

Supplied by :  
Funakoshi Co., Ltd.