

DIRAS2 antibody

Product Information

Catalog No.:	FNab02393
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Displays low GTPase activity and exist predominantly in the GTP-bound form.

Immunogen information

Immunogen:	DIRAS family, GTP-binding RAS-like 2
Synonyms:	GTP-binding protein Di-Ras2 Distinct subgroup of the Ras family member 2 DIRAS2
Observed MW:	22 kDa
Uniprot ID :	Q96HU8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02393(DIRAS2 antibody) at dilution of 1:500