

Product Description: VPAC1 antibody, human, C-term.
VIPR1
Vasoactive intestinal peptide receptor 1
G Protein-coupled receptor
Polyclonal antibody raised against the
C-terminus of the human VPAC1 receptor

BIOTREND Chemikalien GmbH
Im Technologiezentrum Köln
Eupener Str. 157 • D-50933 Köln
Tel. +49 221 949 83 20
Fax. +49 221 949 83 25
jaeger@biotrend.com
www.biotrend.com


Product-No.: S7-9020-S-R

Immunogen: residues 438-457 coupled to KLH (NP_004615)
KLH- TRVSPGARRSSSFQAEVSLV

Host: Rabbit

Form: serum with 0.02% NaN₃ as preservative

Quantity: 50 µl



VPAC1 receptor immunoreactivity in
mucosa of human small intestine.

Applications: ICH (p), IC, W
W = Western blot (1:5000)
IC = Immunocytochemistry (1:5000)
ICH (p) = Immunohistochemistry of formalin-fixed, paraffin-
embedded tissues (1:2500)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4 °C short term
-20 °C long term

