

SLC22A5 antibody

Product Information

Catalog No.:	FNab05966
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	solute carrier family 22(organic cation/carnitine transporter), member 5
Synonyms:	Organic cation/carnitine transporter 2 High-affinity sodium-dependent carnitine cotransporter Solute carrier family 22 member 5 SLC22A5 OCTN2
Observed MW:	75 kDa
Uniprot ID :	O76082

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab05966(SLC22A5 Antibody) at dilution of 1:50



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with FNab05966(SLC22A5 antibody) at dilution of 1:500