

## MAP2K6 antibody

## **Product Information**

Catalog No.: FNab05118

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months(Avoid repeated freeze / thaw cycles.)

## **Background**

This gene encodes a member of the dual specificity protein kinase family, which functions as a mitogen-activated protein(MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases(ERKs), act as an integration point for multiple biochemical signals. This protein phosphorylates and activates p38 MAP kinase in response to inflammatory cytokines or environmental stress. As an essential component of p38 MAP kinase mediated signal transduction pathway, this gene is involved in many cellular processes such as stress induced cell cycle arrest, transcription activation and apoptosis.

## **Immunogen information**

Immunogen: mitogen-activated protein kinase kinase 6

Synonyms: Dual specificity mitogen-activated protein kinase kinase 6 (MAP kinase

kinase 6, MAPKK 6)|MAPK/ERK kinase 6 (MEK 6)|Stress-activated

protein kinase kinase 3 (SAPK kinase 3, SAPKK-3, SAPKK3)|MAP2K6|MEK6|MKK6|PRKMK6|SKK3

Observed MW: 37 kDa Uniprot ID: P52564

**Application** 


Reactivity: Human, Mouse, Rat




Tested Application: ELISA, WB, IHC, FC

Recommended dilution: WB: 1:500-1:2400; IHC: 1:50-1:500

Image:



Immunohistochemistry of paraffin-embedded human breast cancer using FNab05118(MAP2K6 antibody) at dilution of 1:100



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab05118(MAP2K6 antibody) at dilution of 1:500