

MOUSE ANTI ELASTASE 1
MONOCLONAL ANTIBODY

CLONE:	203F4
QUANTITY:	0.1mg
SPECIFICITY:	Human Elastase 1
IMMUNOGEN:	Human Elastase 1
Type/Subtype:	IgG2b κ
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a detection antibody in a human Elastase 1 ELISA combination with the monoclonal capture antibody 207E10.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.

Produced by:
Japan Clinical Laboratories, Inc

Supplied by :
Funakoshi Co., Ltd.