



Product Data Sheet

Product Name: Amylin (1-37), human, amide

Catalog Number: AS-60254-1 (1 mg) **Lot Number:** See label on vial

Sequence: H-Lys-Cys-Asn-Thr-Ala-Thr-Cys-Ala-Thr-Gln-Arg-Leu-Ala-Asn-Phe-Leu-Val-His-Ser-Ser-Asn-Asn-Phe-Gly-Ala-Ile-Leu-Ser-Ser-Thr-Asn-Val-Gly-Ser-Asn-Thr-Tyr-NH₂ (Disulfide bridge: 2-7) (3-letter code)
KCNTATCATQRLANFLVHSSNCFGAILSSTNVGSNTY-NH₂ (Disulfide bridge: 2-7) (1-letter code)

Molecular Weight: 3904.4

Peptide Purity: >95%

Appearance: Lyophilized white powder

Peptide Reconstitution: Amylin (1-37) peptide

Storage: Amylin (1-37) peptide is shipped at ambient temperature. Upon receipt, store lyophilized peptide at -20°C or lower. Reconstituted peptide can be aliquoted and stored at -20°C or lower.

Related Products:

Name	Cat #	Size
Amylin (1-37), human (KCNTATCATQRLANFLVHSSNCFGAILSSTNVGSNTY)	AS-60804	1 mg
Amylin (1-37), rat (KCNTATCATQRLANFLVRSSNNLGPVLPPTNVGSNTY-NH ₂)	AS-60253-05 AS-60253-1	0.5 mg 1 mg
Amylin (8-37), rat (ATQRLANFLVRSSNNLGPVLPPTNVGSNTY-NH ₂)	AS-21749	1 mg

For Research Use Only