

Anti human TR beta 1 mouse monoclonal antibody

TR beta1: Thyroid Hormone Receptor beta1

Code No PP-H3825A-00

Clone No. H3825A

Lot. A-1

Concentration 1 mg/mL

Volume 100 uL

Ig Class G2a

Description Thyroid hormone receptor beta (TRb, THRB; NR1A2) is a member of nuclear hormone receptors. The ligands for TRb are the thyroid hormones, which exist in two forms: T4 (3,5,3',5'-tetraiodo-L-thyronine) and T3 (3,5,3'-tri-iodo-L-thyronine). TRb plays critical roles in the differentiation, growth, metabolism and physiology of a wide variety of tissues. The major partners of TRs are the RXRs, which strongly enhance their ability to bind to specific DNA sequences and contribute to the specificity of action of TR.

Nomenclature NR1A2

Genbank X04707

Origin Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant TR beta 1 (2-100 aa) .

Specificity This antibody specifically recognizes human TR beta 1 but does not recognize TR beta 2 and TR alpha. Not yet tested in other species.

Purification Ammonium sulfate fractionation

Formulation Physiological saline with 0.1% NaN₃ as a preservative.

Application / Recommended Concentration

In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot 1 ug/mL

Non reducing Western Blot Not yet tested

ELISA 0.1 ug/mL

Immunoprecipitation Decide by use

Supershift Assay Not yet tested

Chromatin immunoprecipitation Not yet tested

Immunohistochemistry Not yet tested

Storage Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.

Reference

Notes Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Aug 10, 2006