
TARDBP antibody

Product Information

Catalog No.:	FNab08573
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	TAR DNA binding protein
Synonyms:	TAR DNA-binding protein 43 (TDP-43) TARDBP TDP43
Observed MW:	43 kDa
Uniprot ID :	Q13148

Application

Reactivity:	Human, Mouse, Rat, Zebrafish
Tested Application:	ELISA, IF, WB, IHC, IP, FC
Recommended dilution:	WB: 1:200-1:1000; IP: 1:500-1:1000; IHC: 1:50-1:5000; IF: 1:200-1:800
Image:	



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab08573(TARDBP antibody) at dilution of 1:50



A549 cells(shcontrol and shRNA of TDP43) were subjected to SDS PAGE followed by western blot with FNab08573(TARDBP antibody) at dilution of 1:1000.