



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human STAT1 Polyclonal Antibody, Rabbit

Code # : KB3682GNPAF [FAF01313], (Lot# XXX)

Gene : STAT1 (Signal transducer and activator of transcription 1-alpha/beta, Transcription factor ISGF-3 components p91/p84)

Immunogen : GX5083 (GST-fusion protein, 168 amino acids)
ELQQLDSEKFLQVHQLYDDSFMEIRQYLAQWLEKQDWEHAANDVVSFATIRFHDLLSQLDDQY
SRFSLENNFLLQHNIRKSKRNLDNFQEDPIQMSMIIYSCLKEERKILENAQRFNQAQSGNIQST
VMLDKQKELDSKVRNVKDKVMCIEHEIKSLEDLQDEYDFK

Format : Affinity Purified Rabbit IgG

Constitution : PBS containing with 40% glycerol and 0.02% of NaN₃

Quantity : 50 µg (250 µL)

Specificity : This antibody detects human STAT1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>