

---

## PPL antibody

### Product Information


Catalog No.:	FNab10515
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	periplakin
Synonyms:	Periplakin 190 kDa paraneoplastic pemphigus antigen 195 kDa cornified envelope precursor protein PPL KIAA0568
Observed MW:	204 kDa
Uniprot ID :	O60437

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:50-1:200
Image:	



Immunofluorescent analysis of MCF7 cells using FNab10515(Periplakin antibody) at dilution of 1:50.



Mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab10515(Periplakin antibody) at dilution of 1:1000