



AteloCell[®]

DATA SHEET

Atelocollagen Hanks' Medium

Cat No.: HAN-02

Lot No.: 329031

ORIGIN:	Bovine dermis
CONCENTRATION:	1.9mg/mL
pH:	7.30
STERILITY TEST:	Negative (medium: TGC-I and SCD)
CELL CULTURE TEST:	Normal (Human Fibroblast)
REFERENCES:	Refer to the AteloCell [®] website

FOR RESEARCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.

3-14-3 Mejiro, Toshima-ku, Tokyo Japan TEL: +81-3-3950-6600 FAX: +81-3-3950-6602

<http://www.atelocell.com>