

SOWAHC antibody

Product Information

Catalog No.:	FNab00419
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ankyrin repeat domain 57
Synonyms:	Ankyrin repeat domain-containing protein SOWAHC Ankyrin repeat domain-containing protein 57 Protein sosondowah homolog C SOWAHC ANKRD57 C2orf26
Observed MW:	70 kDa
Uniprot ID :	Q53LP3

Application

Reactivity:	Human
Tested Application:	ELISA, IP, IHC
Recommended dilution:	IP:1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cerebellum tissue slide using FNab00419(ANKRD57 Antibody) at dilution of 1:50



IP Result of anti-ANKRD57 (IP:FNab00419, 3ug; Detection:FNab00419 1:1000) with HeLa cells lysate 3000ug.