

## NeoPower Basic - Basic Electrophoresis Power Supply

#Cat: NB-12-8003



### Features

- ◆ Output: constant voltage, constant current and constant power.
- ◆ Safety: Automatic monitor of over-voltage, voltaic arc, no load, sudden change of load, overload, and short circuit. Leakage protection, function of pause.
- ◆ Others: LCD readable voltage, current, time set. 4-group parallel connection with multi-electrophoresis apparatus.

### Technical parameter

- ◆ Output range : voltage 10-300V ; current: 1-400mA; 75 W
- ◆ Sensitive limit: voltage 1V, current: 1mA
- ◆ Display mode: 192×64LCD
- ◆ Time set: 1min–99h59min
- ◆ Dimension : 300x240x100mm (LxWxH)
- ◆ Weight : 2.0kg