



# AteloCell<sup>®</sup>

## DATA SHEET

**Atelocollagen Bovine dermis 5mg/mL**

**Cat No.: IPC-50**

**Lot No.: 111035**

ORIGIN:	Bovine dermis
CONCENTRATION:	5.2mg/mL
SOLUTION:	1mmol/L HCl
pH:	3.07
STERILITY TEST:	Negative (medium: TGC- I and SCD)
CELL CULTURE TEST:	Normal (Human Fibroblast)
SPECIFIC OPTICAL ROTATION: $[\alpha]_D^{20}$	-430°
PURITY:	Total collagen >95% by SDS-PAGE
REFERENCES:	Refer to the AteloCell <sup>®</sup> website <a href="http://www.atelocell.com">http://www.atelocell.com</a>

FOR RESERCH USE ONLY, NOT FOR HUMAN BODY.

KOKEN Co., Ltd.

3-14-3 Mejiro, Toshima-ku, Tokyo Japan TEL: +81-3-3950-6600 FAX: +81-3-3950-6602

<http://www.atelocell.com>