



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human SMARCC2 Polyclonal Antibody, Rabbit

Code # : KB4171GNPAF [FAF01767], (Lot# XXX)

Gene : SMARCC2 (SWI/SNF complex subunit SMARCC2, SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily C member 2, SWI/SNF complex 170 kDa subunit, BRG1-associated factor 170)

Immunogen : GX5112 (GST-fusion protein, 170 amino acids)
MAVRKKDGGPNVKYYEAADTVTQFDNVRLWLKGKYNKYIQAEPPTNKSLSSLVVQLLQFQEEV
FGKHVSNAPLTKLPIKCFDFKAGGSLCHILAAAYKFKSDQGWRRYDFQNP SRMDRNVEMFMT
IEKSLVQNNCLSRPNIFLCPEIEPKLLGKLDIIRHQGTVED

Format : Affinity Purified Rabbit IgG

Constitution : PBS containing with 40% glycerol and 0.02% of NaN₃

Quantity : 50 µg

Specificity : This antibody detects human SMARCC2 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute

2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>