

IL27RA antibody

Product Information


Catalog No.:	FNab10921
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Interleukin-27 receptor subunit alpha
Synonyms:	Interleukin-27 receptor subunit alpha (IL-27 receptor subunit alpha, IL-27R subunit alpha, IL-27R-alpha, IL-27RA) Cytokine receptor WSX-1 Cytokine receptor-like 1 Type I T-cell cytokine receptor (TCCR) ZcytoR1 IL27RA CRL1 TCCR WSX1
Observed MW:	70 kDa
Uniprot ID :	Q6UWB1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10921(IL27RA Antibody) at dilution of 1:1000