



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Human DR1 Polyclonal Antibody, Rabbit

**Code # :** KB3457GNPAF [FAF01290], (Lot# XXX)

**Gene :** DR1 (Protein Dr1, Down-regulator of transcription 1, TATA-binding protein-associated phosphoprotein, Negative cofactor 2-beta, NC2-beta)

**Immunogen :** GX5089 (GST-fusion protein, 176 amino acids)  
MASSSGNDDDLTIPRAAINKMIKETLPNVRVANDARELVVNCCTEFIHLLISSEANEICNKSEKKTIS  
PEHVIALESLSLGFSGYISEVKEVLQECKTVALKRRKASSRLENLGIPEEELLRQQQLFAKARQQ  
QAELAQQEWLQMQQAAQQAQLAAASASASNQAGSSQDEEDDDDI

**Format :** Affinity Purified Rabbit IgG

**Constitution :** PBS containing with 40% glycerol and 0.02% of  $\text{NaN}_3$

**Quantity :** 50  $\mu\text{g}$  (250  $\mu\text{L}$ )

**Specificity :** This antibody detects human DR1 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at  $-20^\circ\text{C}$ . Avoid freeze-thaw cycles.

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>