



Laboratory pens
for Histology and Cytology

CatNo. E6103KS

Lot:
Storage: Room temperature

Quantity:	12 pens
Application:	The versatile laboratory marker writes on almost all surfaces made of glass, plastics and paper. With this permanent marker, e.g. slides, embedding cassettes, Petri dishes and CDs can be permanently labeled.
Description:	<ul style="list-style-type: none">- ultra-fine writing tip- black writing- writes on slides and embedding cassettes- writes without pretreatment of surfaces- waterproof, odorless and non-toxic- resists most laboratory chemicals- does not dry out and has a long life

E6103KS 230109-1/1

