

ALDOA antibody

Product Information

Catalog No.:	FNab00301
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a key role in glycolysis and gluconeogenesis. In addition, may also function as scaffolding protein(By similarity).

Immunogen information


Immunogen:	aldolase A, fructose-bisphosphate
Synonyms:	Fructose-bisphosphate aldolase A Lung cancer antigen NY-LU-1 Muscle-type aldolase ALDOA ALDA
Observed MW:	39 kDa
Uniprot ID :	P04075

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab00301(ALDOA antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00301(ALDOA antibody) at dilution of 1:500