

COL14A1 antibody

Product Information

Catalog No.:	FNab10495
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis.

Immunogen information


Immunogen:	Collagen alpha-1(XIV) chain
Synonyms:	Collagen alpha-1(XIV) chain Undulin COL14A1 UND
Observed MW:	None
Uniprot ID :	Q05707

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	IHC, IF
Recommended dilution:	IHC: 1:50-1:200; IF: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon carcinoma slide using FNab10495(col14A1 Antibody) at dilution of 1:200



Immunofluorescent analysis of HeLa cells using FNab10495(col14A1 Antibody) at dilution of 1:50.