

A18842

Leader in Biomolecular Solutions for Life Science



# TMBIM6 Rabbit pAb

Catalog No.: A18842

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

27kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

---

Enables endoribonuclease inhibitor activity and ubiquitin protein ligase binding activity. Involved in several processes, including negative regulation of RNA metabolic process; negative regulation of intrinsic apoptotic signaling pathway; and response to L-glutamate. Acts upstream of or within negative regulation of calcium ion transport into cytosol. Located in endoplasmic reticulum membrane and mitochondrial membrane. Biomarker of cervical squamous cell carcinoma and prostate carcinoma.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

7009

### Swiss Prot

P55061

### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMBIM6 (NP\_003208.2).

### Synonyms

BI-1; TEGT; BAXI1; TMBIM6

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.