

MOUSE ANTI PSA
MONOCLONAL ANTIBODY

CLONE:	No.56
QUANTITY:	0.1mg
SPECIFICITY:	Human Prostate specific antigen (PSA) of both free form and complex from.
IMMUNOGEN:	Purified antigen from human seminal fluid
Type/Subtype:	IgG2a κ
APPLICATIONS:	Equimolar –response assay (This antibody can be used as a detection antibody in a human PSA ELISA combination with the monoclonal capture antibody 4D10.)
FORMAT:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.

Produced by:
Japan Clinical Laboratories, Inc

Supplied by :
Funakoshi Co., Ltd.