

MS4A1 antibody

Product Information


Catalog No.:	FNab01441
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	membrane-spanning 4-domains, subfamily A, member 1
Synonyms:	B-lymphocyte antigen CD20 B-lymphocyte surface antigen B1 Bp35 Leukocyte surface antigen Leu-16 Membrane-spanning 4-domains subfamily A member 1 MS4A1 CD20
Observed MW:	33 kDa
Uniprot ID :	P11836

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:2000-1:15000; IHC: 1:1000-1:15000
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01441 (CD20 Antibody) at dilution of 1:10000



Raji cells were subjected to SDS PAGE followed by western blot with FNab01441 (CD20 Antibody) at dilution of 1:8000