


03 Flow Cytometry Antibodies

Advantages

- ▶ High sensitivity/stability, excellent fluorescein labeling technology
- ▶ Multiple fluorescein labeled, hundreds of options
- ▶ More cost-effective, sufficient stock

Multiple fluorescein labeled



Flow Cytometry Colour Matching recommendation

Free Trial

Species	Cell types	Recommendation	Cat.NO
Human	T/B/NKCells	CD45-PerCP	PCP-30039
		CD3-FITC	FITC-30004
		CD16-PE	PE-30061
		CD56-PE	PE-30008
		CD19-APC	APC-30066
Human	Th1/Th2 Cells	CD3-PerCP/Cyanine5.5	PCP55-30004
		CD4-FITC	FITC-30005
		IFN-Y-PE	PE-30053
		IL4-APC	APC-30043
Mouse	Th1/Th2 Cells	CD3-PerCP/Cyanine5.5	PCP55-30002
		CD4-FITC	FITC-30001
		IFN-Y-PE	PE-30074
		IL4-APC	APC-30026
Human	Treg Cells	CD4-FITC	FITC-30005
		CD25-PE	PE-30035
		CD127-APC	APC-30033
Mouse	Treg Cells	CD4-FITC	FITC-30001
		CD25-PE	PE-30017
		FOXP3-APC	APC-30055

New Launch

ELISA Kit/Antibody/Reagent



Wuhan Fine Biotech Co., Ltd.


Distributed by:
CliniSciences Group



Wuhan Fine Biotech Co., Ltd (FineTest®) is an ISO 9001:2015 certified manufacturer which specializes in providing ELISA Kits, Antibodies, Recombinant Proteins and relevant Reagents for research use since 2014.

The factory has a 1200m² GMP workshop covering the core production of main product lines. With professional tech support and rich production experience, our high quality products are widely sold to over 60 countries by cooperating with authorized distributors. Meanwhile we work closely with influential biological platforms like Citeab, Biocompare to promote online marketing.

Now more than 3600+ papers using FineTest's products are cited by international journals and there is an annual increase of about 1000+ citations.



1200m²

10 K class GMP clean room

3600+

Citations in total till now

1000+

New Citations increase annually

60+

Countries across global market

01 QuickTest ELISA Kits



Quick

Time-saving 2-2.5h




Stable

Special protective agent added, higher stability and low CV



Multi-options

2000+ items, various sample validation



Distinctive

Color-coded reagents for avoiding mistakes

Comparison

★ QuickTest ELISA Kits

Time 90-120min

Reagents Color-coded, Ready to use

Sample Size 50ul

VS


Common ELISA Kits

4h

No color-coded and needs to be diluted

100ul

Color change diagram during the experimental process



Selected hot-selling products

Cat.NO	Product Name	Sensitivity	Range
QT-EU3602	ACTH(Corticotropin) ELISA Kit	7.5pg/ml	12.5-800pg/ml
QT-EH0027	Human ACE2 QuickTest ELISA Kit	0.234ng/ml	0.391-25ng/ml
QT-EH2593	Human ADP QuickTest ELISA Kit	0.938ng/ml	1.563-100ng/ml
QT-EH20018	Human CRP QuickTest ELISA Kit	0.938ng/ml	1.563-100ng/ml
QT-EH3058	Human FGF23 QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EH0150	Human GDF15 QuickTest ELISA Kit	14.063pg/ml	23.438-1500pg/ml
QT-EH0410	Human GFAP QuickTest ELISA Kit	0.188ng/ml	0.313-20ng/ml
QT-EH0164	Human IFN γ QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EH3267	Human IL-17A QuickTest ELISA Kit	18.75pg/ml	31.25-2000pg/ml
QT-EH0185	Human IL-1 β QuickTest ELISA Kit	2.344pg/ml	3.906-250pg/ml
QT-EH0201	Human IL-6 QuickTest ELISA Kit	0.28pg/ml	0.46-30pg/ml
QT-EH20020	Human Lipocalin-2 QuickTest ELISA Kit	18.75pg/ml	31.25-2000pg/ml
QT-EH2243	Human MAP1LC3B QuickTest ELISA Kit	0.094ng/ml	0.156-10ng/ml
QT-EH0222	Human MCP-1 QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EH0238	Human MMP-9 QuickTest ELISA Kit	0.188ng/ml	0.313-20ng/ml
QT-EH0018	Human MPO QuickTest ELISA Kit	0.938ng/ml	1.563-100ng/ml
QT-EH1205	Human NEFL QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EH3468	Human OC/BGP QuickTest ELISA Kit	0.141ng/ml	0.234-15ng/ml
QT-EH4007	Human PEA QuickTest ELISA Kit	18.75pg/ml	31.25-2000pg/ml
QT-EH0022	Human PLGF QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EH3832	Human sTREM-1 QuickTest ELISA Kit	18.75pg/ml	31.25-2000pg/ml
QT-EH0287	Human TGF β 1 QuickTest ELISA Kit	18.75pg/ml	31.25-2000pg/ml
QT-EMK0189	Monkey TSH QuickTest ELISA Kit	0.0094uIU/ml	0.0156-1uIU/ml
QT-EM0109	Mouse IL-1 β QuickTest ELISA Kit	1.875pg/ml	3.125-200pg/ml
QT-EM0112	Mouse IL-2 QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EM0121	Mouse IL-6 QuickTest ELISA Kit	0.938pg/ml	1.563-100pg/ml
QT-EM1592	Mouse IL-8/cxcl15 QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EM0153	Mouse PCSK9 QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-EM0205	Mouse VEGF QuickTest ELISA Kit	9.375pg/ml	15.625-1000pg/ml
QT-ERB0126	Rabbit TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	4.688pg/ml	7.813-500pg/ml



02 Cell Function Assay Kits

FineTest® offers a comprehensive of various cell function assay kits, including detection for cell apoptosis, cell cycle, cell proliferation/cytotoxicity/viability etc. We provide a complete detection solution for cell function and help you to conduct cell-related researches.



Free Trial

Our self-developed cell function assay kits are featured with high sensitivity, short assay time, high accuracy and wide choices. Welcome scientific researchers and institutions to contact your local distributor to order.

Cell apoptosis is different from cell necrosis in cell morphology and physiology. Changes on cell membrane, cytoplasm, nucleus, nucleic acid and apoptosis-associated proteins are obvious. Compared with normal or necrotic cells, characteristics are also distinct.

FineTest® can offer various detection methods for cell proliferation/cytotoxicity/viability, including EdU, Calcein AM/PI, CCK-8 assay etc. You can choose the relevant cell assay kit based on the sample type and detection purpose.

Using the property that intracellular DNA can bind to fluorescent dyes (e.g. propidium iodide-PI), the DNA content of cells varies from period to period and thus binds to different fluorescent dyes, and the fluorescence intensity detected by the flow cytometer varies.


Enhanced Cell Counting Kit 8 (CCK8/WST-8)

MTT-Cell Based Proliferation /Toxicity Assay Kit

Cell Viability/Cytotoxicity Kits

EdU Proliferation Kits

Cell Cycle Kits



▲ Obvious cells grouping . Half reagent doses can get good results.

▲ Obvious cells grouping and accuracy result

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: osterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: danmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

